

## Product datasheet for RN212796

### Pfn2 (NM\_030873) Rat Untagged Clone

#### Product data:

Product Type: Expression Plasmids  
 Product Name: Pfn2 (NM\_030873) Rat Untagged Clone  
 Tag: Tag Free  
 Symbol: Pfn2  
 Mammalian Cell Selection: Neomycin  
 Vector: pCMV6-Entry (PS100001)  
 E. coli Selection: Kanamycin (25 ug/mL)  
 Fully Sequenced ORF: >RN212796 representing NM\_030873  
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGGCCGGTTGGCAGAGCTACGTGGATAACCTGATGTGCGATGGCTGCTGCCAGGAGGCCGCCATTGTCG  
 GCTACTGCGACGCCAAATACGCTGCGGCAGCCACGGCCGGGGCGTCTCCAGAGCATCACGCCAGCAGA  
 AATAGATGTGATTATAGGAAAAGACCGGGAAGGTTTCTTTACCAATGGTTTGACTCTTGGAGGAAAGAAG  
 TGCTCTGTGATCAGAGATAGCCTATACGTTGACTCGGACTGCACAATGGACATCCGGACAAAGAGTCAAG  
 GTGGGGGCCAACATAACAACGTTGCTGTTGGCAGAGCTGGGAGAGTCTTGGTCTTTGTAATGGGAAAAGA  
 AGGGGTCCATGGAGGCGGATTGAATAAGAAGGCATACTCAATGGCAAAATACTTGAGAGACTCTTGGGTT  
 CTAGCTGCTAGGCAGACTGTTAAGTATTAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI  
 ACCN: NM\_030873  
 Insert Size: 450 bp  
 OTI Disclaimer:

Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_030873.1</a></u> , <u><a href="#">NP_110500.1</a></u>
<b>RefSeq Size:</b>	1966 bp
<b>RefSeq ORF:</b>	450 bp
<b>Locus ID:</b>	81531
<b>UniProt ID:</b>	<u><a href="#">Q9EPC6</a></u>
<b>Cytogenetics:</b>	2q26
<b>Gene Summary:</b>	binds actin; regulates the rate of actin polymerization [RGD, Feb 2006]