

## Product datasheet for **RN212801**

### **Cnot7 (NM\_001107313) Rat Untagged Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Cnot7 (NM\_001107313) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Cnot7  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN212801 representing NM\_001107313  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGCCAGCAGCAACTGTAGATCATAGCCAAAGAATTTGTGAAGTTGGGCTTGTAACCTGGATGAAGAGA  
TGAAGAAAATCCGTCAAGTTATCCGAAAATAATTATGTTGCTATGGACACCGAGTTCCAGGCGTTGT  
TGCAAGGCCCATTTGGAGAATTCAGAAGCAATGCTGACTATCAGTACCAACTGTTGCGATGTAATGTAGAC  
TTGTTAAAGATAATTCAGCTTGGACTGACATTTATGAATGAGCAAGGAGAATACCCTCCAGGAACATCAA  
CTTGCCAGTTTAAATTTAAATTTAATTTGACGGAGGACATGTATGCTCAGGACTCTATAGAGCTACTGAC  
AACGTCTGGCATCCAGTTTAAAAACATGAGGAGGAAGGGATTGAGACCCAGTATTTGCAGAACTTCTT  
ATGACCTCAGGAGTGGTCTCTGTGAAGGGTCAAATGGTTATCATTTTCATAGTGGTTATGACTTTGGCT  
ATTTAATCAAAATTCGACCAACTCTAACTTGCCTGAAGAAGAGCTTGATTTCTTTGAGATCCTTCGATT  
ATCTTCCCTGTCATCTATGATGTGAAGTACCTCATGAAGAGCTGCAAGAACCTCAAGGGTGGACTGCAG  
GAAGTTGCTGAGCAGTTAGAGCTGGAGCGGATAGGACCTCAACACCAGGCAGGATCTGACTACTGCTTA  
CAGGAATGGCCTTTTCAAATGAGAGAAATGTTCTTTGAAGATCACATTGATGATGCCAAATACTGTGG  
CCACTTATACGGCCTTGTTCTGGCTCATCTATGTACAGAACGGCACAGGGAATGCATATGAAGAGGAA  
GCCAGCAAGCAGTCG**TGA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001107313  
**Insert Size:** 858 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_001107313.1, NP_001100783.1</u>
<b>RefSeq Size:</b>	2614 bp
<b>RefSeq ORF:</b>	858 bp
<b>Locus ID:</b>	306492
<b>Cytogenetics:</b>	16q12.1