

Product datasheet for RN212839

Cxcr4 (NM_022205) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
 Product Name: Cxcr4 (NM_022205) Rat Untagged Clone
 Tag: Tag Free
 Symbol: Cxcr4
 Synonyms: MGC108696
 Mammalian Cell Selection: Neomycin
 Vector: pCMV6-Entry (PS100001)
 E. coli Selection: Kanamycin (25 ug/mL)
 Fully Sequenced ORF: >RN212839 representing NM_022205
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGAAATATACACTTCGGATACTACTCCGAAGAAGTAGGGTCTGGAGACTATGACTCCAACAAGGAAC
 CCTGCTCCGGGATGAAAACGAAAACCTCAACAGGATCTTCTGCCACCACATCTATTTTATCATCTTCTT
 GACTGGCATAGTGGGCAATGGTGGTAATCCTGGTCAATGGGTTACCAGAAGAAGCTGAGGAGCATGACA
 GACAAGTACCGGCTGCACCTGTCCGTGGCTGACCTCCTCTTTGTCATCACACTCCCCTTCTGGGCAGTGG
 ACGCCATGGCTGACTGGTACTTTGGGAAATTTTATGTAAGGCTGTGCATATCATCTACACCGTCAACCT
 TTACAGCAGTGTCTCATCCTGGCCTTCATCAGCCTGGACCGCTACCTTGCCATTGTCCACGCCACCAAC
 AGCCAGAGGCCGAGGAAGCTGCTGGCTGAAAAGGCCGTCTATGTGGTGTCTGGATCCCCGCCCTCCTCC
 TGACTATCCCTGACATCATCTTCGCGGATGTCAGCCAGGGGACGGCAGGTACATCTGTGACCGCCTTTA
 CCCCAGACGCTGTGGATGGTGGTGTCCAGTTCCAGCACATCATGGTGGTCTCATCCTGCCGGGCATC
 GTCATCCTGTCTGTTACTGCATCATCTCCAAGCTGTCACTCCAAGGGCCACCAGAAGCGCAAGG
 CCCTCAAGACTACGGTCACTCTTATCCTGGCTTTCTTGCCTGCTGGCTACCGTATTACGTGGGGATCAG
 CATCGATTCTTCATCCTTTGGAGTCAAGCAAGGATGTGAGTTCGAGAGCGTCGTGCACAAGTGG
 ATCTCCATCAGGAGGCCCTCGCCTTCTCCACTGTTGCCTGAACCCATCCTACGCCCTCCTCGGGG
 CCAAATTCAGAGCTCCGCGCAGCATGCACTCAATTCCATGAGCAGAGGCTCCAGCCTCAAGATCCTTTC
 CAAAGGAAACGGGGTGGACTCTTCCGTCTCCACAGAGTCAGAATCCTCAAGTTTTCACTCCAGCTAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Chromatograms: https://cdn.origene.com/chromatograms/ja2689_e04.zip



[View online »](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_022205
Insert Size:	1050 bp
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_022205.3 , NP_071541.2
RefSeq Size:	1726 bp
RefSeq ORF:	1050 bp
Locus ID:	60628
UniProt ID:	O08565
Cytogenetics:	13q13
Gene Summary:	G protein-coupled receptor for the chemokine stromal cell-derived factor-1; may have a role in brain development [RGD, Feb 2006]