

Product datasheet for **RN212888**

Clec2g (NM_001048075) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids

Product Name: Clec2g (NM_001048075) Rat Untagged Clone

Tag: Tag Free

Symbol: Clec2g

Synonyms: CLEC2D11; Clr-b; Clr11; Clrb; RGD1563148

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >RN212888 representing NM_001048075
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGAGTGCCAAAAAGGCTTCCCAACCTATGCTCAATACCACAGGCAGCCTGCAGGAGGGAGAAATGGGTA
 AAATGTTCCATGGGAAATGTCTCAGAATCGTCTCCCCGAATCTCCTGCTAAGCTTTACTGCTGCTATGG
 AGTGATCATGGTCCTCAGTGTAGCTGTAGTTGCTCTTTCTGTTGCTTTGTCAGTAAAAATGACACCACAG
 ATCTCAACCATAAATACCTACGCTGCTTGCCCGAGAACTGGATTGGAGTTGAAATAAATGTTTTATT
 TTTCTGAATATGCAAGTAACTGGACATTCAGCCAGACCTTCTGCAAGGCACAAGAGGCCGAGCTAGCACG
 ATTTGACACCGAGGAGGAGCTGAATTTCTTAAGCAGATACAAAGGGAGTTTGTATTACTGGATCGGCCTG
 CACAGAGAGTCCTCAGAGCACCTTGGAAGTGGACAGACAACACTCAGTATAACTACTCGCTTTCATCC
 GGGGAGTGGAAGATATGCCTACCTGAACGACATCGGGATCAGCAGTGCCAGGGTCTATGCAGACAAAAG
 ATGGAGCTGTAGCAAACTTAACAGCTATAGCCTCCAATGCAAACTCCTTTTCTCCTATG**TAG**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_001048075

Insert Size: 624 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001048075.1, NP_001041540.1</u>
RefSeq Size:	624 bp
RefSeq ORF:	624 bp
Locus ID:	362447
UniProt ID:	<u>Q0H8B9</u>
Cytogenetics:	4q42
Gene Summary:	Receptor for KLRB1B that protects target cells against natural killer cell-mediated lysis. [UniProtKB/Swiss-Prot Function]