

## Product datasheet for **RN213008**

### **Cbx6 (NM\_001012119) Rat Untagged Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Cbx6 (NM\_001012119) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Cbx6  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN213008 representing NM\_001012119  
**Red**=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGAGCTGTCTGCAGTAGGCGAGTGGGTCTTCGCGGCCGAATCCATCATCAAACGCCGGATCCGAAAGG  
GACGCATCGAGTACCTGGTGAATGAAAGGGTGGCGATCAAGTACAGCACTTGGGAGCCAGAGGAGAA  
CATTCTGGATTCGAGGCTCATCGCAGCCTTCGAGCAGAAGGAAAGGAAACGTGAGCTGTATGGGCCAAAG  
AAGAGGGGACCTAAACCCAAAACCTTTCCTCCTAAAGCCAAGCGCCAGCGCCTCCTCGCCCAAGCTGCACT  
CCAGCGTGCAGTGCACCGGCTCAAGAAAGACATCCGCCGCTGCCACCGTATGTCCCGCCGCCCCCTGCC  
ACGACCAGACCCACAAGGGGTAACCTGGCCTGCGCCACCCATCTCGCCCTTTCTGAGACCGTGCGC  
ATCATCAACCGCAAGGTGAAGCCCGGGAACCAAGCGGAACCGCATCATCCTGAACCTGAAGGTGATTG  
ACAAGGGCCCGGAGGTGGCAGCACTGCCAGGGCAGCGGGCACTGGCCCGCCGAAAGTCCCCTCGC  
GAACCGGGTCAATGGGAAGAGCAAGAAGTTCAGCGAGAGCATGCTCCGACCCAGATCCGCCACATGAAG  
TTTGGCACCTTTGCACTCTACAAGCCCCCTCCCGCCCTCTGGTGCCTCTACTGCGGGCAAAGCTGACG  
TGGCCTCCTCGGGCCTGGGCTGCTCCTGGCCACCCAGCTGCTGCCCTTTGACGCACACAGCTCTAG  
CTCTCCGGTGCCTCACCACACTGCACTCCTGACCCAGATGATGCCACCCAAAGCTCCTCCCC  
GAGACCATAAGCCGATCTGCGTCCAACCTGGCGAGAGTCAAGAGTCTCGATCTGTCCATCCCTCCTGAGG  
CAGCTGCCACCGCCAGCGGGTGCCTCCTGATGCTCACTGCTGCTGCGGGCCAGGCCTACACACAGTTCT  
GGAGGCCACAGGTACCGGCTCCTCCGAGCCTGAGGCTGGGGACTGGCGCCCGAGATGTCGCCATGCTCC  
AATGTAGTCGTACAGATGTCAACAGCAACCTCCTGACCGTGACCATCAAGGAATTCTGCAGCCCTGAGG  
ATTCGAGAAGGTGGTCTGCTGGGGTGGCAGGTGCTACAGGGGGTGGGGTGGCACTGGGCCAAGTAAGTG  
A

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI



[View online »](#)

<b>ACCN:</b>	NM_001012119
<b>Insert Size:</b>	1191 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001012119.1</a></u> , <u><a href="#">NP_001012119.1</a></u>
<b>RefSeq Size:</b>	3095 bp
<b>RefSeq ORF:</b>	1191 bp
<b>Locus ID:</b>	315136
<b>Cytogenetics:</b>	7q34