

## Product datasheet for RN213140

### Sh3gl1 (NM\_031239) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Sh3gl1 (NM\_031239) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Sh3gl1  
**Synonyms:** SH3P8  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN213140 representing NM\_031239  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGTCGGTGGCGGGCTGAAGAAGCAGTTCTACAAGCGAGCCAGCTGGTCAGCGAGAAAGTTGGTGGG  
CCGAAGGGACAACTGGATGATGACTTCAGAGAGATGGAAAAGAAAGTGGATATCACCAGTAAGGCCGT  
GGCAGAGGTGCTGGTCAGAACCATAGAATATCTGCAACCTAACCCAGCCTCGAGGGCCAAGTCTGCTGG  
CTGAATACTGTATCCAAGATCCGGGGCCAAGTGAAGAATCCTGGTACCCACAGTCAGAGGGTCTGCTGG  
GAGAGTGCATGGTCCGCCATGGCAAGAACTAGGCGGAGAGTCCAACCTTTGGCGATGCTCTGCTAGATGC  
AGGTGAGTCTATGAAGCGCCTGGCTGAGGTGAAGGACTCACTGGACATCGAGGTCAAGCAGAACCTCATC  
GACCCACTGCAGAACCTGTGTGACAAGGATCTGAAGGAGATCCAGCACCCTGAAGAAGTTGGAGGGCC  
GCCGCTTGACTTTGACTACAAGAAGAAGCGCCAGGGCAAGATCCCTGATGAGGAGCTGCCTCAGGCCCT  
AGAGAAGTTTGAGGAGTCCAAGGAGGTGGCGGAGACCAGTATGCACAACCTCCTGGAGACTGATATTGAG  
CAGGTGAGCCAGCTCTCAGCCCTAGTGGATGCCAGCTGGACTACCACCGCAGGCAGTGCAGATCCTGG  
AGGAGCTGGCTGATAAGCTGAAGCGCAGGGTGAAGAAGCTCCTCAGCCCCAGGAGGGAGTTCAAGCC  
CAGGCCCCAGGAGCCCTTTGAGCTTGAGAGCTGGAGAGCTGGAGCAGCCCAATGGGGGATTCCCTGTGCCTCAGCA  
CCCAAGATCACAGCTTCGTCATATTTAGATCAGGGGACAAGCCACCAGGACGCCAGCAAGAGTATGC  
CACCCCTGGACCAGCCAAGCTGCAAGGCCTGTATGACTTTGAGCCAGAGAACGATGGCGAGCTGGGCTT  
CCGAGAGGGTGACCTCATCAGCTTACCAACCAGATCGATGAGAAGTGGTATGAGGGGATGCTGCATGGC  
CAGTCAGGCTTCTCCCACTTAGCTATGTGCAGGTGCTGGTGCCTCTGCCTCAGTGA

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_031239



[View online »](#)

|                               |   |
|-------------------------------|---|
| <b>Insert Size:</b>           | 1107 bp   |
| <b>OTI Disclaimer:</b>        | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).  |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).  |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol> |
| <b>RefSeq:</b>                | <u>NM_031239.2, NP_112518.1</u>   |
| <b>RefSeq Size:</b>           | 2010 bp   |
| <b>RefSeq ORF:</b>            | 1107 bp   |
| <b>Locus ID:</b>              | 81922   |
| <b>UniProt ID:</b>            | <u>O35964</u>   |
| <b>Cytogenetics:</b>          | 9q12  |
| <b>Gene Summary:</b>          | interacts with synaptojanin [RGD, Feb 2006]   |