

## Product datasheet for RN213150

### Blvrb (NM\_001106236) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Blvrb (NM\_001106236) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Blvrb  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Fully Sequenced ORF:** >RN213150 representing NM\_001106236  
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGGCTGTCAAGAAGATCGCGATCTTCGGTGCCACCGAAGGACCGGGCTCACCACACTGGCGCAGGCGG  
 TGCAAGCAGGTTACGAGGTGACGGTGCTGGTTCGAGATTCCAGCAGGTTACCGTCAGAGGGGCCCCAGCC  
 AGCCCACGTGGTGGTGGGAGATGTTCTGCAGGCGGGCGATGTGGACAAGACTGTGGCTGGACAGGATGCC  
 GTCATCGTGCTACTGGGCACTGGCAACGACCTCAGTCTACTACAGTAATGTCCGAGGGCACCCGGAACA  
 TTGTGGCAGCCATGAAGGCACATGGAGTGGACAAGGTCGTGGCCTGCACCTCGGCCTTCTACTCTGGGA  
 CCCGAGCAAGGTGCCCCCTCGCCTACAGGACGTGACCGATGACCACATCCGGATGCATAAGATTCTGGAG  
 GAGTCAGGACTGAAGTACGTGGCTGTAATGCCCCACACATTGGGGACCAACCACTAACCGGGCCCTACA  
 CAGTGACCCTGGATGGACGAGGGCCCTCGAGGGTCATATCCAAGCATGATCTGGGTCACCTCATGTACTAG  
 GTGCCTCACCACCATGAGTATGACGGACAGAAGACCTACCCTTCCCACCAGTATGACT**AG**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** Sgfl-Mlul  
**ACCN:** NM\_001106236  
**Insert Size:** 621 bp  
**OTI Disclaimer:**

Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_001106236.1, NP_001099706.1</u>
<b>RefSeq Size:</b>	884 bp
<b>RefSeq ORF:</b>	621 bp
<b>Locus ID:</b>	292737
<b>Cytogenetics:</b>	1q21