

## Product datasheet for RN213157

### Hbegf (NM\_012945) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Hbegf (NM\_012945) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Hbegf  
**Synonyms:** Dtr; GFHB; Hb-egf; Hegfl  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN213157 representing NM\_012945  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGAAGCTGCTGCCGTCGGTGGTGCTGAAGCTCTTTCTGGCCGAGTGTTGTCCGCGTTGGTGACCGGTG  
AGAGTCTGGAGCGGCTTCGGAGAGGTCTGGCGGCAGCAACCAGCAACCCTGACCCTCCACTGGAACCAC  
AAACCAGCTGCTACCCACGGGAGCTGATCGCGCTCAGGAGTCCAGGACTTGAAGGGACCGATCTGGAC  
CTTTTCAAAGTTGCTTTCTCCTCCAAGCCACAAGCCCTGGCCACCCAGGAAAAGAAAAGAACGGGAAAA  
AGAAGAGGAAAGGCAAGGGGTTAGGAAAGAAGAGATCCATGCCTTAAGAAATACAAGGACTACTGCAT  
CCACGGAGAGTGAGATACCTGAAGGAGCTCCGTATTCCCTCGTGCCACTGCCTCCCTGGTTACCATGGA  
CAGAGGTGTCATGGGCTGACCCTACCGGTAGAGAACCCCTGTACACATATGACCACACTACCGTCTTGG  
CTGTGGTGGCTGTAGTACTGTCTGTCTTCTTGTTCATCGTGGGACTTCTCATGTTACAGGTACCA  
TAGGCGAGGAGGTTATGACTTGGAAAGTGAGGAGAAAGTGAAGTTGGGCATGGCTAGCTCCCACTGA

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_012945  
**Insert Size:** 627 bp  
**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_012945.1</a></u> , <u><a href="#">NP_037077.1</a></u>
<b>RefSeq Size:</b>	1550 bp
<b>RefSeq ORF:</b>	627 bp
<b>Locus ID:</b>	25433
<b>UniProt ID:</b>	<u><a href="#">Q06175</a></u>
<b>Cytogenetics:</b>	18p11
<b>Gene Summary:</b>	activates epidermal growth factor receptor mediated signaling; contributes to neuronal survival; may play a role in neurogenesis during ventral midbrain development [RGD, Feb 2006]