

Product datasheet for **RN213325**

Setsip (NM_194353) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Setsip (NM_194353) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Setsip
Synonyms:	RGD727855; Set; Setl1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN213325 representing NM_194353 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGGCCCCGAAGCGGCAATCTGCGATCCCGCTCAGCCCAAGAAACCCAGACCCATTGCTGCCCCGAAGT
TGGAGGACAAGTTGGCCTCTCCAGCCTGCCGAAGGGAGAAAAAGAACAGCAAGAACTATTGACCATAT
TGATGAAGTACAAAATGAAATGGACAGACTTAATGAACAACCCAGTGAGGAAATTTGAAAGTAGAACAA
AAATATAACAGACTCTGCCAACCATTTTTTTTTTTCAGAAGAGGTCAGAATTGATCGCCAAATCCCAATT
TTTGGGTAAACAATTTGTCAACCATCCACAAGTGTCTGCACCGCGTGGGAGGAAGACGAGGAGGCCCT
GCATTATTTGACCAGAGTTGAAGGGACAGAATTTGGAGACATTAATCAGGTTACAGAATAGATTTTTAT
TTTGATGAGAATCCTTACTTTGAAAATAAAGTTCTCTCAAAGAATTTCTCTGAATGAGAGCGGTGACC
CGTCTTCAAAGTCCACCGAAATCAAATGGAAATCTGGGAAGGATTTGACAAAATGCTCAAGTCAAACGCC
GAATAAGGCCAGCAGGAAGAGGCAGCGTGAAAAGCCAGAGAGTTCTTTCCCTGGTTTACTGACCGTTCT
GACGCAGGTGCTGACGAGTTAGGAGAGGTCGTCAAGGATAATTTGGCCAAACCCCTTGCGGTCTTATT
TGGTTCGGATATGGGTGATGAAGGAGAGGCAGAGGATGATGATGATGACGACGACGAAGAAGAACCGTT
GGAAGATATTGATGAAGGAGACAGGGATGAAGGTGAAGAAGATGACGATGAGGAGGAAGGGGAGGAAGGA
GGACGAAGCGAGGATGA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_194353
Insert Size:	858 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_194353.1, NP_919334.1</u>
RefSeq Size:	858 bp
RefSeq ORF:	858 bp
Locus ID:	290432
Cytogenetics:	15q21
Gene Summary:	protein phosphatase regulator protein; human homolog appears to be involved in histone acetylation [RGD, Feb 2006]