

Product datasheet for RN213344

Spaca3 (NM_001105820) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Spaca3 (NM_001105820) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Spaca3
Synonyms:	Lyzc
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN213344 representing NM_001105820 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAAGCTAGGAGCCGGGCTCCAGAAGACAGCCGTGCCCGCTGGGATCACTTGGCTGGCCCTGGCCT
ATCTGCTCAGCTGCCTGCTTGCCTCCAGCAAGGCCAAGGTCTTCAGTCGCTGTGAGCTGGCCAAAGTGCT
GCATGACTTCGGTCTGGAAGGCTACAGGGGATATAACCTGGCTGACTGGATCTGCCTTGCCTACTACACA
AGCGGCTTCAATACGGATGCTGTGGATCACGAAGCTGACGGAAGCACCAACAACGGCATCTTCCAGATCA
GCAGCCGGAAGTGGTGC AAAACCTCGCCCGAACGGCCCAATCTTTGCAGGATATACTGCACGGATTT
GTTGAGCAATGATCTCAAAGATTCTGTCGCTGTGT CATGAAGATAGCGCAAGAGCCCAGGGCCTGGGC
TACTGGGAATCCTGGAAGCACCATTGTCAGGGCAGGGACCTCAGTGACTGGTGGATGGCTGTGACTTCT
AG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001105820
Insert Size:	492 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001105820.1, NP_001099290.1</u>
RefSeq Size:	939 bp
RefSeq ORF:	492 bp
Locus ID:	287557
Cytogenetics:	10q26