

## Product datasheet for RN213348

### Ubfd1 (NM\_001034911) Rat Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Ubfd1 (NM_001034911) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Ubfd1
Synonyms:	RGD1306614; Ubph
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN213348 representing NM_001034911 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCGCGGCCGGAGCCCCGGATGGCATGGAGGAGCCAAGCATGGACACGGAGGCCGAGGCTGTGGCCA  
CCGAGGCACCCGACGGCCCCCTCAACTGCGTGGAGGCCGAAGCAGTCGTAGGGACAGCGCGGAGGACTC  
CTGCGACGCGCGAGGCAACCTGCAGCCGGCCCCGGCCAGCCCCCTGGGGACCCTGCGGCCAGGCCTCG  
GTCAGCAACGGTGAGGACGCGGGTGGCGGTGCGGGCAGAGAGCTGGTGGATCTGAAGATAATCTGGAACA  
AGACCAAGCACGACGTGAAGTTCCTCTGGACAGCACGGTTCGAGCTAAAGCAGAAGATTCACCTCCAT  
TACAGGTCTCCCACCCGCCATGCAGAAAGTCATGTATAAGGGACTTGTTCTGAAGACAAGACGTTGAGA  
GAAATAAAAGTGACCAAGTGGGGCCAAGATCATGGTGGTTGGCTCCACTATAAATGATGTTTTAGCAGTAA  
ACACCCAAAAGATGCTGCCAGCAGGATGCAAAGGCTGAAGAAAGCAAGAAGGAGCCCCCTCGCAGGCA  
GAAGCAACACAGGAAAGTGTGGACAAAGGAAAACCTGAAGATGTGATGCCATCTGTTAAGGGTGCCAG  
GAGCGCTGCCAACGGTACCCTTGTCTGGCATGTACAATAAGTCTGGAGGGAAAGTGAGACTGACCTTTA  
AGCTAGAGCAAGACCAGCTGTGGATCGGCACTAAAGAGCGGACTGAGAACTGCCATGGGCTCCATTAA  
AAATGTTGTGAGTGAACCTATTGAAGGACATGAAGACTACCACATGATGGCGTTTCAGTTGGGCCCCACG  
GAAGCCTCTTACTACTGGGTGTACTGGGTTCCAACCTCAATATGTGGATGCAATCAAAGACACTGTGCTGG  
GCAAGTGCCAGTATTTTGA

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001034911



[View online »](#)

<b>Insert Size:</b>	930 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001034911.1</a></u> , <u><a href="#">NP_001030083.1</a></u>
<b>RefSeq Size:</b>	1235 bp
<b>RefSeq ORF:</b>	930 bp
<b>Locus ID:</b>	293454
<b>Cytogenetics:</b>	1q36