

Product datasheet for RN213399

Scgb3a1 (NM_001013180) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Scgb3a1 (NM_001013180) Rat Untagged Clone
Tag: Tag Free
Symbol: Scgb3a1
Synonyms: Hin1; PnSP-2; Pnsp2; secretoglobin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN213399 representing NM_001013180
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGAAGCTCACCGCCATCTTTCTAGTGTTCTGTGTGGCTCTGCTCAGCGACTCTGGCGTTGCTTTCTTCG
 TGGATTCAGTGGCAAGCCTGTGGTAGAACCCGTGGCTGCCCTTGCTCCAGCTGCAGAGGCTGTGGCAGG
 GGCTGTGCTAGCCTACCATTAAAGCCACTTGCCATCCTGAGGTTTCATCCTGGCCAGCCTGGGCATCCCA
 TTGGATCCTCTCATAGATGGTTCAGGAAGTGCCTCACCGAGCTGGGCCCTGAGGCTGTAGGAGCTGTGA
 AGTCACTGCTGGGGGCCCTGACAACGTTTCGGTTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-MluI

ACCN: NM_001013180

Insert Size: 315 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



View online »

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001013180.1](#), [NP_001013198.1](#)

RefSeq Size: 331 bp

RefSeq ORF: 315 bp

Locus ID: 338418

Cytogenetics: 10q21

Gene Summary: secretoglobin; may have a role in airway function [RGD, Feb 2006]