

Product datasheet for RN213430

Exog (NM_001106866) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Exog (NM_001106866) Rat Untagged Clone
Tag: Tag Free
Symbol: Exog
Synonyms: Endogl1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN213430 representing NM_001106866
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGCTGCAAAGGGCTTCGCTTCGCGTCTCAGGGACTCGCGGAGGTTTCTGAATGGCTTCTTGCCGGGG
CTGTAGTGGCGCCGAGGAGCAGGGCTGACGGCGCTGCAGTTCCTCCAGCGTCGGGACCGGGGTCGC
GAGGCCGCGAGGCAGCCGACGAATCTGTAGAAAAAGTTGTCTTGAACAGTTTGGATTCTTTAACT
GGAACAGAGACCAGGCGTTACACTAACCATGCCTTGTCTGTATGATCAGGCAAAGCGGGTCCAGATGGG
TTCTTGAACATATATCCAAGGACAAGATCATCGGTGACGCGGACAGAAAGCATTGTAATTAAGCCGGA
TCCTACTGTCCCTCAGCCTTCAGCGCCCTCAATGAAGACTATATTGGAAGTGGCTGGTACAGGCCAC
ATGGCCCCAGCAGGAAATAACAAGTTTTCCAGTGAAGCCATGGCAGAAACCTTTTACCTTTCAACATCG
TGCCTCAGAAATTTGAGAATAATTCCGGCTACTGGAATAGGATAGAAATGTAAGCTGAGCTGACGGA
GAGATTTGAAGATGTTTGGATAGTATCGGGACCTTAACTTACCTCATACTAGAAACGATGGAACGAAG
ACAGTTAGTTACCAGGTGATCGGTGAGGACAATGTGGCAGTCCCCTCACACCTTTACAAGGTAATCCTGG
CTCGCAGAAGCCAGAATCTACCGAACCACTGGCACTGGGGCCCTTTGTGGTACCAAATAAAGCCATCGG
GTTCCAGTCCCAGTTAAGCGAGTTCAGGTGAGCCTCCATGACCTGGAGAAGATGTCAGGCCTGGTGTTT
TTCCCCCATCTGGATAGAACAGGGACATCCGGAACATTTGCTCTGTGGACACCTGCAAACCTCTGGGGT
TCCAGGAGTTCACTCTGTACCTGAGCACAAGGAAGATCGACGGAGCCCGTCTGTGGCCAGACTGGAGAA
GGTTCTGGAAGCTTTGAAGAGCTCTGGGGTGAACCTGACGATTACTTCTGAGTCGCTATGGGAAAAAG
CTGGAGGAACTGAAGGCTAAGGAGCAGAAGGACGCCAGCTGGAAAAGCAGTCC**TAG**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAAGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001106866



Insert Size:	1107 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001106866.1, NP_001100336.1</u>
RefSeq Size:	4085 bp
RefSeq ORF:	1107 bp
Locus ID:	301062
Cytogenetics:	8q32