

Product datasheet for **RN213474**

Pih1d3 (NM_001024890) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Pih1d3 (NM_001024890) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Pih1d3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN213474 representing NM_001024890 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAATCAGAAAATGCAGAAGCTGAGAACATAGAAACAGGGAAGTTGGAAACAAAGACGTGGAAGCTA
AAAATATGCGATTTGAGACACTGTCTTCTGCCTCGGCCCTGCAAGCACTGTCCAATCTCCTGAATCCTGA
AGAAGAGGATGATTTTGATTTTCAACAGGGACAATGCTCATCTACTATTGGAGCCATGGGCCCTGGGAAT
ATTGGACCACCTAAAGCAAAAGAGTCCAAAGCCATTCTGAACCCAGGAGTGATGAAAGTGAGAACATCT
GGAATCCAGAAGAGGTTCCAGAAGGAGCAGAGCATGATGATATATGGGATGCTCGAGAAGTCCCAGATTA
TGAGATTGTGTTCCAACAAACGGTGGAAGTGAAGACGTGTACCTAGGATTAACAAGAAAGGATCCTTCA
ACAGCATGCTGCCAGGAGCTGGTGGTGAATAAACTGCCGATACAAACCCTTCTGAAATTCAAATTG
ACATCCAAGAAATGCTCCTTGATCTTCGTACTCTCTAAAGAGAGCTGCTAGTGAACCTTCCCAGCCTGT
GGAACATAACAGTGCTAAAGCGTCCTATATCTGGGAAGCTGAAACACTTGAAGTCAGGATGACCCTCCAA
AGGATTTGGATTTTAATATTTC**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001024890
Insert Size:	657 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001024890.1, NP_001020061.1</u>
RefSeq Size:	1600 bp
RefSeq ORF:	657 bp
Locus ID:	315915
Cytogenetics:	Xq32