

## Product datasheet for **RN213489**

### Stk32b (NM\_001107224) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Stk32b (NM\_001107224) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Stk32b  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN213489 representing NM\_001107224  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGGCGGGAACCACTCCCACAAGCCCCAGTGTTTGATGAGAACGAAGAAGTCAACTTTGACCACTTCC  
AGATCCTTCGGGCGATTGGTAAAGGGAGTTTCGGAAAGGTATGCATCGTGCAGAAGCGAGACACAAAGAA  
GATGTATGCCATGAAGTACATGAACAAGCAGAAGTGCCTGGAGAGGGATGAGGTGCCAACGTGTTCCGG  
GAGCTGCAGATCATGCAGGGATTAGAACACCCCTTCTGTGTGAACCTGTGGTACTCCTCCAGGATGAGG  
AAGACATGTTTATGGTGGTGGACCTGCTGCTGGGCGGGGACCTGCGCTACCACCTACAGCAGAACGTGCA  
CTTCACAGAGGGAGCCGTGAAGCTCTACGTCTGTGAGCTGGCCCTGGCACTGGAGTACCTGCAGAGGTAC  
CACATCATCCACAGAGACATCAAACCAGACAACATCCTGCTGGATGAGCATGGGCATGTGCACATCACGG  
ACTTCAACATCGCCACAGTCTGAAAGGAAGTGAAGAGGCTTCTCCATGGCTGGCACAAGCCCTACAT  
GGCTCCTGAAATCTTTCAGGTGTATGTGGATGGAGGCCCTGGTACTCATACCCCGTGGACTGGTGGTCC  
CTGGGTGTACAGCCTATGAGCTGCTACGAGGCTGGAGGCCATATGAAATCCACTCGGCCACACCCGTCG  
ACGAGATCCTCAACATGTTCAAGGTGGAACGTGTCCACTACTCCTCCACGTGGTGTGAGGGGATGGTGGAA  
CCTGCTGAAGAAGCTGTGACCAAGGACCCGAGATCCGTCTGTCCAGCCTCCGTGACATCCAGAGCATG  
ACCTACCTGGCTGACATGAACTGGGATGCAAGTGTCAATAAGGCGCTGATGCCTGGCTTTGTGCCAATA  
AGGGGAGACTGAACTGTGACCCACATTTGAACCTGAAGAGATGATTCTAGAATCCAAGCCGCTTCACAA  
AAAGAAGAAGAGGTTGGCCAAGCATAGGTCCAGAGACAGCACGAAGGACAGCTGTCCCCTGAATGGACAC  
TTGCAGCAGTGTGGGAGACAGTGCAGGAGGAGTTCATCATATTAACCGAGAGAAAACCTCAGGAGGCAGC  
AGGGACATGATGGTCACTCTCAGACCTGGATGGCCGAGTAGGAAACCAGACCTCGAGCAAGGTGCAAGA  
CGGTCGAAACAACAACATCCTGACCCATACATGCCCTCGAGGCTGCAGCAGCTAA

**ACGCGT**ACGCGGCCGCTCGAGCAGAAAACCTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI



[View online »](#)

<b>ACCN:</b>	NM_001107224
<b>Insert Size:</b>	1245 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001107224.2</a></u> , <u><a href="#">NP_001100694.1</a></u>
<b>RefSeq Size:</b>	1245 bp
<b>RefSeq ORF:</b>	1245 bp
<b>Locus ID:</b>	305431
<b>Cytogenetics:</b>	14q21