

Product datasheet for RN213502

Prmt6 (NM_001106466) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Prmt6 (NM_001106466) Rat Untagged Clone
Tag: Tag Free
Symbol: Prmt6
Synonyms: Hrmt116
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN213502 representing NM_001106466
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGC**C

ATGTCGCTGAGCAAGAAAAGAAAGCTTGAGTCGGGGTTCGGCGGCAGGAGGGAGGGAGCTGAGGAGG
 AAAATGGCGGGGAGCAGGAGGCAGCCCCGCCACGACCCCGGAGGACCAAGCGCGAGCGGGACCAGCTGTA
 CTACGAGTGCTACTCCGACGTATCGGTCCACGAGGAGATGATCGCCGACCGCGTCCGACCGACGCCTAC
 CGCCTAGGCATCCTGAGGAAGTGGGCCGCTCTGCGAGGCAAGACCGTGTGGACGTGGGCGCAGGCACCG
 GCATTCTTAGCATCTTCTGCGCCAGGCCGGGGCCCGCGCTATACGCTGTGGAGGCCAGCGCCATCTG
 GCAACAGGCCAGGAGGTGGTGGCGCTCAACGGGTTGGAGGACCGGTGCACATCCTGCCAGGCCGGTG
 GAGACGGTGGAGCTGCCGGAGCAAGTGGACGCCATCGTAAGCGAGTGGATGGGCTACGGACTACTGCACG
 AGTCCATGCTGAGCTCCGTGCTCCACGCGCGGACCAATGGTGAAGGAGGGCGGTCTCCTCCTGCCAGA
 TTCCGCTGAGCTCTTCTGGCCCCGATTAGCGACCAGATGCTGGAATGGCGCTTGGGTTTCTGGAGCCAA
 GTGAAGCAGCACTATGGCGTGGATATGAGCTGCATGAAAAGTTCGCCACGCGCTGCCTCATGGCCATT
 CGGAGATCGTGGTGCAGGGTCTGTCCGGAGAGGACGTGCTGGCCCGGCCACAGCGCTTTGCCAGCTCGA
 GCTGGCCCCGCGCTGGCCTGGAGCAGGAGCTGGAGCCGGTGTGGCGGGCGCTTCCGCTGCAGCTGCTAT
 GGCTCCGCGCCTTACATGGTTTTGCCATCTGGTTTTAGGTAACCTTTCTGGAGGGGACTCGGAGAAAC
 CTCTGGTGTGTCCACCTCACCTTTCCACCCGGCCACGCACTGGAAGCAGGCGCTCCTCTACTTGAATGA
 GCCGGTCCGGTGAACAAGATACAGACATTTCTGGAGAGATCACCTGCTGCCCTCCCGGGACAACCC
 CGGCGTCTGCGCTACTTCTGCGCTACAAAGTGGGGACCATGAGGAAAAGACCAAGACTTTGCCATGG
 AGGACTGA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



[View online >](#)

ACCN:	NM_001106466
Insert Size:	1128 bp
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).</p>
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001106466.1, NP_001099936.1</u>
RefSeq Size:	3312 bp
RefSeq ORF:	1128 bp
Locus ID:	295384
Cytogenetics:	2q41