

Product datasheet for **RN213534**

Cryz11 (NM_001013044) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Cryz11 (NM_001013044) Rat Untagged Clone
Tag: Tag Free
Symbol: Cryz11
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN213534 representing NM_001013044
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGAAAGGCTTGTATTTCAACAGAGTCTACAAATGAAGAAGTAACATTTGTATTTCAAGAAAAGGAAA
 ATGTCTCTGTACAGAGGATAAATTTGTGAGAGTGCAAGTTAAAGCCTGTGCCCTGAGCCACATAAACAC
 AAAGCTTCTGGCAGAGATGAAGATGAAAAGGAATTCCTCCCTGTTGGGAGAGAAGTCGCTGGGATTGTA
 TTAGATGTTGGCAAGAAGGTATCGTTCCTTCAGCCAGATGATGAAGTCGTCGGGATTTTGCCTTTGGATT
 CGGAAGACCCCGACTCTGCGAAGTTATCCGAGTACAGGAGCATTACCTGGTTCATAAGCCAGAAAAGGT
 GCCATGGACGGAAGCCGCTGGAGTCATCCGGGATGGAGTCCGAGCCTACACAGCCCTCTATTACCTTTCT
 CAGCTCTCTCCCGGAAGTCCGTGCTGATCATGGATGGGGCAAGTGCATTTGGCACAATAGCCATTCAGT
 TAGCACACCACAGAGGAGCCAAGGTGATCTCAACAGCACACAGCCTGGAGGACAAGCAACATCTTGAAAG
 GCTTAGGCCTTCTAGCCCGAGTGATTGATGTGTCTAGTGGGAAAGTCCATGTTGCTGAAAGGTGCTTG
 GAAGAAACAGGTGGCCTGGGAGTAGATATTGTCATAGATGCAGGAGTGAGATTATATAGTAAAGATGATG
 AACCAGCTGTAAAACAACAACACTACCACATAAGCATGACATCATCACACTTCTCAGTGTGGAGGCCATTG
 GATAACGACAGAAGAAAACCTTCAGTTGGATCCTCCAGATAGCCATTGCCTTTTCCCTCAAGGGAGCAACA
 GTGGCTTTCCCTAATGATGAAGTTTGGAAATTTGTCAAATGCACAACAGGGAAAATATCTTTGCATCTTAA
 AAGACGTGATGGAGAAGCTATCAGCTGGTGTGTTGAGACCTCTGCTTGTGATGAACCCATTCCACTGTATGA
 GGCCAAAGTTTCCATGGAAGTTGTTGAGAAAAGTCAAGAAAGAAAAAAGCAAGTTGTTCAATTT**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001013044
Insert Size: 1047 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001013044.1</u> , <u>NP_001013062.1</u>
RefSeq Size:	1853 bp
RefSeq ORF:	1047 bp
Locus ID:	288256
Cytogenetics:	11q11