

Product datasheet for **RN213546**

Zcrb1 (NM_001034940) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Zcrb1 (NM_001034940) Rat Untagged Clone
Tag: Tag Free
Symbol: Zcrb1
Synonyms: RGD1309851
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN213546 representing NM_001034940
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGAGCGGTGGCTTGGCTCCAAGTAAGAGCACAGTTTATGTATCCAACCTGCCTTTCTCCCTGACAAACA
ATGATTTATACCGGATATTTTCCAAGTATGGCAAAGTTGTAAGGTTACTATAATGAAAGACAAGGATAC
TAGGAAGAGTAAAGGTGTTGCATTCATTTATTTTGGATAAGGACTCTGCCTTAACTGCACCAGAGCA
ATAAATAATAACAGTTATTCGGTAGAGTGATAAAGGCAAGCATCGCCATTGACAACGGAAGAGCAGCCG
AGTTCATCCGGAGGCGAAACTACTTTGATAAATCTAAGTGCTATGAATGTGGGAAAAGTGACACTTAAG
TTATGCCTGTCCTAAAAACATGCTTGGAGAACGTGAGCCTCAAAGAAGAAAAGAAAAGAAAAGAAAAG
AAAGTCCCTGAACCGGAAGAAGAAATTGAGGAAGTAGAAGTGAGTGAAGAAGAAGGGGAGGACCCGGCTC
TCGACAGTCTCAGTCAGGCCATAGCATTCCAGCAGGCCAAAATTGAAGAAGAACAAAACAAATGGAGACC
CAGCACAGGGGGTCCCTCGACCTCCGATGATTCTAGAAGACCAAGAATAAAGAAGAGCGCATACTTCAGT
GATGAGGAGGAACCTCAGTGACT**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001034940
Insert Size: 654 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001034940.1</u> , <u>NP_001030112.1</u>
RefSeq Size:	1326 bp
RefSeq ORF:	654 bp
Locus ID:	362990
UniProt ID:	<u>Q499V6</u>
Cytogenetics:	7q35