

Product datasheet for **RN213647**

Limd2 (NM_001025715) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Limd2 (NM_001025715) Rat Untagged Clone
Tag: Tag Free
Symbol: Limd2
Synonyms: RGD1309967
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >RN213647 representing NM_001025715
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGTTCCAGGCTGCCGGAGCCGCCAGGCCACCCCTCTCATGAAGCCAAAGGCAGCAGTGGCAACAGCA
CGGTGCAGCGGTCTAAGTCCTTTAGCTTGC GGCTCAGGTGAAGGAGACCTGTGCAGCCTGCCAGAAGAC
TGTGTACCCTATGGAGCGGCTGGTGGCAGACAAGCTCATTTTCCACAACCTTTGTTTCTGTTGCAAACAC
TGCCACACCAAACCTCAGCCTGGGCAGTTACGCTGCAATGCACGGTGAATTTACTGCAAACCTCACTTTC
AGCAGCTGTTTAAGAGCAAAGGCAACTACGATGAGGGGTTTGGTCGTAACAGCACAAGGAGCTCTGGGC
CCACAAGGAGGTGGACTCAGGCACCAAAACAGCCT**G**A

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-MluI
ACCN: NM_001025715
Insert Size: 387 bp
OTI Disclaimer:

Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001025715.1</u> , <u>NP_001020886.1</u>
RefSeq Size:	1187 bp
RefSeq ORF:	387 bp
Locus ID:	360646
UniProt ID:	<u>Q4KM31</u>
Cytogenetics:	10q32.1
Gene Summary:	Acts as an activator of the protein-kinase ILK, thereby regulating cell motility. [UniProtKB/Swiss-Prot Function]