

Product datasheet for RN213678

II15 (NM_013129) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: II15 (NM_013129) Rat Untagged Clone
Tag: Tag Free
Symbol: II15
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN213678 representing NM_013129
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGAAATTTTGAACCATATATGAGGAATACATCCATCTTGACTACCTGTGTTTCCTTCTCAACAGTC
ACTTCTTAAGTGGCTGGCATCCATGTCTTCATTTGGGCTGTGTCAGTGTAGGCTCCCTAAAACAGA
GGCCAAGTGGATAGATGTAAGATACGATCTGGAGAAAATTGAAAGTCTTATTCAATTTATCATATTGAT
ACTACTCTATACACTGACAGTGACTTTTCATCCCAGTTGCAAAGTTACTGCAATGAAGTCTTCTCCTGG
AATTACAGGTTATTTGCACGAGTACAGTAACATGACTCTTAATGAAACAGTAAGAAACGTGCTCTACCT
TGCAAACAGCACTCTGTCTTCTAACAAGAATGTAATAGAGTCTGGCTGCAAGGAATGTGAGGAGCTGGAG
GAGAGAACTTCACGGAGTTTTTGCAGAGTTTTATACATATTGTCCAATGTTTCAACACGTCTGA

ACGGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_013129
Insert Size: 489 bp
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_013129.2](#), [NP_037261.1](#)

RefSeq Size: 489 bp

RefSeq ORF: 489 bp

Locus ID: 25670

UniProt ID: [P97604](#)

Cytogenetics: 19q11

Gene Summary: human homolog is a cytokine involved in the regulation of T and natural killer cell activation and proliferation [RGD, Feb 2006]