

Product datasheet for **RN213994**

Taar2 (NM_001008512) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Taar2 (NM_001008512) Rat Untagged Clone
Tag: Tag Free
Symbol: Taar2
Synonyms: Gpr58
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN213994 representing NM_001008512
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGACATCTTTGAAGCCCAACAGGAAACCTTCGATTGTTTCAAGATATGGGAATGGATCTTGCCAGAGA
ATGAAAGGTCTCTGGGGTCCGAGCTGCAATGTACTCACTGATGGCAGGAGCCATCTTCATCAGCATCTT
TGGCAACCTGGTCATGATCATTTCCATTTCTACTTCAAGCAGCTTCACACACCGACCAACCTCCTCATT
CTCTCCATGGCTGTCACCGACTTCTGTTGGGATTCCACATCATGCCATACAGTATGGTCCGGTCAGTGG
AAAAGTCTGGTATTTTGGGCTTACCTTTTGAAGATCCATTATAGCTTTGACCTGATGCTTAGCATAAC
GTCCATTTTCCACCTTGGCTCGGTGGCCATAGATAGATTTTACGCCATCTGTCACCTTTGCATTATTGC
ACCAAAATGACTATTCGGTCTGCAAGCGGCTGCTGCTCGTCTGCTGGTCAAGTCCCGGGGGCTTTGCCT
TTGGGGTGGTCTTTTCAAGGCTTATGCTGATGGGATTGAAGGCTATGACATTTTGGTTGCATGTTCCAG
TTCCTGCCAGTCATGTTTAAACAAGCTCTGGGGGACCACCTTATTTGTGGCGGGCTTTTACTCCGAGC
TCCATGATGGTGGGATTTATGGCAAAATTTTCGAGTATCCAAAAACATGCTCGCGTAATTGACAAC
TGCCAGAAAAATCAAATAATCAAATGAGGAAGGCAAAAAGGCAGCAAAAACCTTGGGGATAGTGATGGG
TGTTTTCTACTGTGTTGGTTTCCATGTTTCTTACGATCCTGCTAGATCGTTTTCTGAATTTCTCTACG
CCTGCAATTCTGTTGATGCCTAACTTGGTTCGGCTATTTTAATTCCACATGTAATCCCTTAATTTATG
GTTTCTTCTACCCCTGGTTTCGACAGGACTCAGGTACATTTTGTGGTAAAATTTTCAGCTCACATTT
CCACAACACTAATTTGTTTACGCAAAAAGAAACAGAA**TAG**

ACGGTACGGCGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001008512



[View online »](#)

Insert Size:	1020 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001008512.1</u> , <u>NP_001008512.1</u>
RefSeq Size:	1020 bp
RefSeq ORF:	1020 bp
Locus ID:	294121
UniProt ID:	<u>Q5QD25</u>
Cytogenetics:	1p12
Gene Summary:	Orphan receptor.[UniProtKB/Swiss-Prot Function]