

Product datasheet for **RN214118**

Fbxl4 (NM_001107919) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Fbxl4 (NM_001107919) Rat Untagged Clone
Tag: Tag Free
Symbol: Fbxl4
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN214118 representing NM_001107919
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCACCAGTGTTCATGTTAACAGTCTGACCATGTTTTATTATATGTGCCTTCGGCGTCGAGCCA
GGACAGCTACGAGAGGAGATATGATGAGCAGCCATAGGACTATAGTATCAAACAGCCGGACTTCCCCTCT
CAATGCAGAGGTGGTCCAGTATGCTAAGGAAGTGGTGGATTTTCAGTTCCTACTACGGCAGCGAGAACAGC
ATGTCCTATACCATGTGGAATTTGGCAGGTGTACCAAATGCTTCCCAAGTCTGGCGACTTCACGCAGA
CGGCTGTGTTTCGAACCTACGGAACATGGTGGGATCAGTGTCTAGTGCATCCTTGCCATTGAGAAGAAC
GCCGTCAGTTTTGAGGCCAGGACTATGTGGAAGTACTTTGAGCAACAAGTGTATCCCACAGCTGTC
CATGTGCTGGAGACCTACCATCCTGGAGCAGTCATCAGAATTCTGGCCTGCTCGGCAAAATCCCTACTCCC
CAAATCCTCCAGCAGAAGTAAGATGGGAGACTCTTGGTCCGAGAGACCCACGAAGGTGAATGCTTCCCA
GGCTCGCCAGTTAAACCTCGTATTAAGCAAATAAATTTCCCACAAACCTTATACGACTGGAAGTAAAC
AGTTCTCTCTGGATTATTACACCGAATTAGATGCAGTTGTGCTACATGGTACGAAGGACAAGCCACTGC
TCTCTCAAGACTGCACTTGTAGACATGAATGATCTAGAAGATGATGACTATGAAGAGAAGGACGGCTG
TGAAATGGACGCTCTTAACAAAAAGTTTAGCAGTGTACCTTGGGGATGGGCCAAGTAACGGGTATTTT
GATAAACTACCATATGAGACCTGTGTAGATT**AG**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-MluI
ACCN: NM_001107919
Insert Size: 873 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001107919.1, NP_001101389.1</u>
RefSeq Size:	2480 bp
RefSeq ORF:	873 bp
Locus ID:	313101
Cytogenetics:	5q21