

Product datasheet for RN214146

PrI5a2 (NM_001044702) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids

Product Name: PrI5a2 (NM_001044702) Rat Untagged Clone

Tag: Tag Free

Symbol: PrI5a2

Synonyms: Ghd2; Prlpp

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >RN214146 representing NM_001044702
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCAGCTGTCTTTGACTCGACCACACTCCTGGGCACTCCTGATACTGCCATTGTCAAACCTGCTTCTGT
 GGGAGAATGTGGTCTCTTTGCCCATGTGTGAAGTAATAAATGGGCAATGCCAGCTGACTCTTGAGAACTT
 GCTTAATCAGGCACACAGCATATCTGAGAACCAATCATCTCACTTCTGAAATATTCAATGAGTTTGTAT
 AAAGAATACTCTTCAGTTTCAGGCTTCTCAGACACAATCCCCCTGGTGTGCTCCAATTATCCCTCCCTG
 TTTCAAATATCACAATACAAGAACAAGAGATTCAGTCTGAGGTCCTCTTGAAAGTGACCATAGATACACT
 CCTTGGCTGGAACAACACCCTGAGTGAGGTAGCAGCAGATATGGATGATCTGAAATCAATTCCTGGTGTT
 GGAGCTTTCATATCTAATGTGAGAAAAATAGCTGCAAAATTCACAAGACTCACAAGTGACTAAAAGAAG
 TCAAGTCATTGCTCAAATTGGTTCGTCTTGAACCTTGAGGAAGATGAGGACAACCCCTGCCTCAGCTGGACT
 TCCATCATTGCATCTGCTAAATAATCCCTCTCGCCTTACCTTCTATCACGCCCTTCTGGGCTGCGTGAGT
 TATGATGCAGAAAAGATTGCAAGTCATGTGAAGTCTGAGGTGCAAAATGATCCATGAAAAATGCT**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_001044702

Insert Size: 699 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001044702.1, NP_001038167.1</u>
RefSeq Size:	919 bp
RefSeq ORF:	699 bp
Locus ID:	691469
Cytogenetics:	17p12