

Product datasheet for **RN214169**

Fam188a (NM_001106122) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Fam188a (NM_001106122) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Fam188a
Synonyms:	RGD1309605
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >RN214169 representing NM_001106122
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**

ATGTCGAAGTGACCAAAGAGCTGCTGGAGCTGGTGTGGGACCAAGAGCAGTCCCGACTGTCGGACA
 CCATCTTCTGCCGCTGGACGCAAGGGTTTGTGTTTAGTGACTCAGAGGGATCTGCATTAGAGCAGTTTGA
 AGGTGGCCCTGTGCTGTTATTGCACCTGTTCAAGCATTCTTTTGAAGAAGCTTCTGTTTTCTCCGAG
 AAGTCGTCCTGGCGGATTGCTCAGAGGACGAGCAGAAGGAGCTCCTTTGTCATACCTTGTGTGATATT
 TAGAAAGTGCTTGTGACAGCTCTGGATCATACTGCATGGTTTCATGGTTGAGAGGAAGGACGACTGAGGA
 CGCCGCTAGGATTTCTGGGAGTCTGCACAGTCTAGTTGCCAAGTGGGACATTCTCTGCCTTGGCTGTC
 GAAGAGCTTGGCTTTGAGCGATTTATGCATTAATTCAAAAAAGATCATTAGAAGTGTATCAGAATTAA
 AAGATGCTGTCTGGACCAGTATTCAATGTGGGAAACAAATTTGGAGTCTTGCTTTTCTCTATTCTGT
 ATTACTGACAAAGGCATTGAAACATAAAGAATTCAATTGAAGATGCAAACGAGCCTTTGATAGACCCT
 GTGTACGGGCATGGCAGCCAAAGCTTAATAAACCTCCTCCTGACGGGACGCGGTCTTAACGTGTGGG
 ATGGCGATAGAGAGTGCTCAGGAATGCAACTTCTCGGTATCCATGAACAAGCAACTGTAGGATTCTTAAC
 ACTAATGGAAGCTTTACGCTATTGTAAGGTTGGTCTTACTTGAAATCTCCAAATTCCTATTTGGATT
 GTTGGCAGCGAACTCACCTCACCGTATTTTTTGCCAAGGATATGGCTTTAGTTGCCCTGAAGCGCCTT
 CGGAACAAGCCAGAAGAGTTTTTCAGACCTATGACCCAGAAGATAATGGATTATAGCCGACACACTTCT
 AGAGGATGTATGAAAGCCTTGGACCTTGTTCAGATCCCGAGTATATAAATCTCATGAAGAACAATTG
 GATCCAGAAGGATTAGGAATCATATTGTTGGGTCCATTCTTCAAGAATTTTTCTGACCAGGGCTCCA
 GTGGTCCAGAATCTTTACTGTCTACCACTACAATGGATTGAAGCAGTCAAATTATAATGAAAGGTAAT
 GTATGTGGAAGGACTGCAGTGGTTATGGGTTTTGAAGACCCCATGCTACAGACAGACGACTCCTATT
 AAACGCTGTCTGCAAAACCAATGGCCGTACATTGAGTTACTGTGGACCACAGACCCTCTCTTCACTCA
ACTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_001106122

Insert Size: 1335 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_001106122.1, NP_001099592.1
RefSeq Size: 2321 bp
RefSeq ORF: 1335 bp
Locus ID: 291320
Cytogenetics: 17q12.3