

Product datasheet for **RN214212**

Cops7b (NM_001108807) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Cops7b (NM_001108807) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Cops7b
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN214212 representing NM_001108807 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCAGCAGTGGGAGGCTTCGGGGGAACAGCACCAGTGACGGGTGCGCGGAGACAGAAAAGTGCCGGAC
GCCGGTCCGATCATGGAAGCTTGACAACCTGCAGGCAGGCTCTGGGAGGCCGAGACATCGGCGAAGAGAA
CAGAGTCGGCGGGGACAGATCTCAAGACCAGAGAATGGCAGGTGAACAGAAACCTCAAGTAACCTCTTG
GAGCAGTTTATTTATTAGCCAAAGGCACCAGCGGCTCAGCCCTCACCCTCTCATAAGCCAGGTACTGG
AGGCTCCTGGAGTTTATGTTTTGGAGAACTACTGGAGTTGGCCAATGTTTCAGGAGCTTGCAGAAGGAGC
TAATGCTGCATATTTGCAGTTGCTGAACCTGTTTGCCTATGGTACATACCCGATTACATAGCCAACAAG
GAGAGCCTGCCAGAACTGAGCGCAGCCAGCAGAACAAAGCTGAAGCATCTTACCATCGTGAGTTTGGCAT
CAAGAATGAAGTGTATCCCCTACTCTGTGCTGCTGAAGGACCTGGAGATGAGGAATCTCCGGAACTAGA
AGACCTTATCATTGAGGCTGTCTACACTGACATCATCCAGGGCAAGCTGGACCAGCGAAACCAGCTCCTG
GAAGTGGATTTCTGCATTGGCCGTGACATCCAGAAGAAAGATATCAATAATATTGTCAAGACCTTGCATG
AATGGTGTGACGGCTGTGAGGCAGTTCTCCTGGGCATTGAGCAGCAAGTTCTGAGAGCCAACCAAGTACAA
AGAGAACCACCACCGACTCAGCAGCAGGTGGAAGCCGAGTTAGCAACATTAAGAAGACACTCAAAGCA
ACAGCGTCCTCCTCAGCTCAGGAAATGGAGCAGCAACTGGCTGAACGGGAGTGTCCCCCATACTGAGC
AGAGGCAACCCACTAAGAAGATGTCAAAGTGAAGGTCTGGTCTCCAGCCGCCAC**TAG**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001108807
Insert Size:	969 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001108807.1</u> , <u>NP_001102277.1</u>
RefSeq Size:	2496 bp
RefSeq ORF:	969 bp
Locus ID:	363273
Cytogenetics:	9q35