

## Product datasheet for RN214215

### Nhp2l1 (NM\_212515) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Nhp2l1 (NM\_212515) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Nhp2l1  
**Synonyms:** Otk27; Snu13  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN214215 representing NM\_212515  
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**GCGATCGCC**

ATGACGGAGGCTGATGTGAACCCAAAGGCCTATCCCCTCGCAGATGCCACCTACCAAGAAGCTGCTGG  
 ACCTTGTTTCAGCAGTCGTGTAACAAGCAGCTTCGAAAGGAGCCAATGAAGCCACCAAACCTCAA  
 CAGAGGCATCTCTGAGTTCATTGTGATGGCAGCAGACGCTGAGCCCTGGAGATCATCCTGCACCTCCCT  
 CTGCTGTGTGAAGACAAGAATGTCCCCTATGTATTTGTACGCTCCAAGCAGGCTTTAGGACGGCCCTGTG  
 GAGTCTCCAGGCCGGTCATCGCCTGTTCTGTAACCATCAAAGAAGGCTCTCAACTAAAGCAGCAGATCCA  
 ATCCATTCAGCAGTCTATTGAGCGGCTCCTGGT**G**TAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** Sgfl-MluI  
**ACCN:** NM\_212515  
**Insert Size:** 387 bp  
**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).  
**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_212515.2](#), [NP\\_997680.1](#)

**RefSeq Size:** 632 bp

**RefSeq ORF:** 387 bp

**Locus ID:** 300092

**UniProt ID:** [P55770](#)

**Cytogenetics:** 7q34

**Gene Summary:** Involved in pre-mRNA splicing as component of the spliceosome. Binds to the 5'-stem-loop of U4 snRNA and thereby contributes to spliceosome assembly. The protein undergoes a conformational change upon RNA-binding.[UniProtKB/Swiss-Prot Function]