

Product datasheet for RN214293

Sirt6 (NM_001031649) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Sirt6 (NM_001031649) Rat Untagged Clone
Tag: Tag Free
Symbol: Sirt6
Synonyms: MGC114368
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN214293 representing NM_001031649
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGTCGGTGAATTATGCAGCAGGGTTGTCGCCATACGCGGATAAGGGCAAGTGCGGGCTGCCGAGATAT
TCGACCCACCAGAGGAGCTGGAATGCAAGGTGTGGGAGCTGGCGGGCTGATGTGGCAGTCTCCACTGT
GGTTTTCCACACCGGCGCAGGCATCAGCACCGCCTCCGGCATCCCCGACTTCAGGGGCCCCATGGCGTG
TGGACCATGGAGGAACGAGGCCTCGCCCCAAGTTTGACATCACCTTTGAGAATGCCCGGCCCTCGAAGA
CCCACATGGCCCTTGTTAGCTGGAACGCATGGGCTTCCTCAGCTTCCTGGTCAGCCAGAATGTAGACGG
GCTACATGTGCGCTCGGGCTTCCCCAGGGACAAGCTGGCCGAGCTGCACGGAAACATGTTTGTAGAGGAG
TGCCCCAAGTGTAAAGACGCAGTACGTGAGAGACACGGTCGTGGGTACCATGGGCCCAAGGCCACTGGCC
GGCTCTGCACCGTGGCCAAGGCGAGGGGACTTCGGGCCTGTAGAGGGGAGCTGAGAGACACCATTCTGGA
CTGGGAGGATGCGTTGCCTGACCGGGACCTAACGCTCGCTGATGAGGCCAGCAGGACTGCAGACCTGTCT
GTCACCTGGGTACCTCACTGCAGATCCGCCCCAGCGGGAACCTGCCCTTGCCACTAAGCGCCGCGGAG
GCCGCTGTGGTCATTGTCAACCTGCAACCCACAAAACATGACCGCCAGGCTGATCTGTGCATCCATGGCTA
CGTGGATGAGGTGATGTGCAAGCTCATGAAGCACCTGGGGCTGGAGATCCCCACCTGGGATGGACCCCGA
GTGCTAGAGAAAAGCCCTGCCACCGCTGCCGCGCCAGTAGCACCCAAGGCTGAGCCCCCTGTGCATCTCA
ATGGTTCTATAAGCCCAAGCCCGACAGCCCTGTACCCACAGGCCTCCCAAAAGAGTGAAGACAGAGGC
TGCCGCCAGCTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001031649



[View online »](#)

Insert Size:	993 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001031649.1</u> , <u>NP_001026819.1</u>
RefSeq Size:	1656 bp
RefSeq ORF:	993 bp
Locus ID:	299638
Cytogenetics:	7q11