

Product datasheet for **RN214382**

118bp (NM_053374) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: 118bp (NM_053374) Rat Untagged Clone
Tag: Tag Free
Symbol: 118bp
Synonyms: lgifbp
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN214382 representing NM_053374
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGACCATGAGACTGTGGCTGTGCAGCAGACCCAGTCCTTGGTGGGTCTGCTTCTATATGTCCATG
TTGTCATTCTGGCCAGAGCCACATCTGCACCTCTGACAACAGCCACTGTCTTAACCAAGAAGCTCAAAGGA
CCCGTGCTCTTCTGGTCTCCAGCAGTCCCAACTAAGCAGTACCCAACATTGGATGTGATCTGGCCAGAG
AAAGAAGTGCCACTGAATGGAACGCTGACCTTGTCTGCACTGCCTGCAGCCGCTTCCCCAACTTCAGCA
TCCTCTACTGGCTGGCAATGGCTCCTTATTGAGCACCTCCCAGGCCGGCTGAGGGAGGGCCACACAAG
CCGGGAGCAGAGGAACGCAAGCACCTGGCTACACAGGGCCTTGGTGTGGAGGAAGTGGCCCTCCCTA
CTAAGTACCAACTTCTCTGTCTGTTTGTGGATCCCGGACAGGTGGCCAGTATCATGTCTTCTGGCCC
AGCTCTGGGATGGATTGAAGACAGCCCCATCCCCTTCTCAAGAAACCTCTCCAGCCACAGCCCAGGACC
CAGATCAGCAGGCCAGGGGCCCAT**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-MluI
ACCN: NM_053374
Insert Size: 588 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_053374.2</u> , <u>NP_445826.2</u>
RefSeq Size:	951 bp
RefSeq ORF:	588 bp
Locus ID:	84388
Cytogenetics:	1q32
Gene Summary:	antagonist for the proinflammatory responses induced by endogenous IL-18 [RGD, Feb 2006]