

Product datasheet for RN214454

Olr1278 (NM_173129) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
 Product Name: Olr1278 (NM_173129) Rat Untagged Clone
 Tag: Tag Free
 Symbol: Olr1278
 Mammalian Cell Selection: Neomycin
 Vector: pCMV6-Entry (PS100001)
 E. coli Selection: Kanamycin (25 ug/mL)
 Fully Sequenced ORF: >RN214454 representing NM_173129
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCAACAGGAAACTATTGCGTATTTCTGAGTTCATTCTGACTGGGCTCTCCAAAAATCAGAACTTC
 AGATGCCTCTCTTTGTCCTCTTCTAGGAATATACATAGTCACTGTAGTGGGGAATTTGGGCATGATTAC
 ACTGATTAGGCTCAGTTCCTCCTGCACACACCAATGACTATTTCTCAGCAGTTTGCCTTTATTGAT
 CTCTGCCATTCCACAGTCATTACCCCAAAATGTTGGTGAACCTTGTGGCAGAGAAGAACATTATCTCT
 ACCTGGATGCATGACTCAGCTTTTTCTTTCTCATTTTTGTATTGCAGAGTGCATATGTTAGCTGC
 AATGGCATATGATCGCTATGTTGCCATCTGTAACCCACTGCTTTACAATGCAATCATGTCCTATCAGAGT
 TACATTTCCATGATCTCTGGAGTCTATATTATTGGTGTGGTTTGCATCAGCTCACACAGGCTTCATGA
 TAAGAAGTCAGTTCTGTAACCTTAGATGTGATCAACCACTATTTCTGTGATCTTCTCCCCTGCTGGAGCT
 TGCACACTCTAGTACTTATGTTAATGAATTGTAATTTTAACTGTGGGACATGTAACATTGTTGTTCCA
 ACCTTGACTATTCTACTTCTTATCTTCATTATTGCCACCATCTACACATTCGTTCTACTGAGGGCA
 GGTACAAAGCATTTCAGTACCTGCAGCTCCACATCTTGGCTGTTGCTGTCTTCTTTGGGTCTGCTGCATT
 TATGACTTACAGCCATCATCAGTCAGCTCTATGGACCAAGGAAAGTGTCTCTGTGTTTTATACCATC
 GTTGTGCCATGCTGAACCCCTTGATCTACAGTCTGAGGAATAAGGATGTTAGTACTGCCCTGAAGAAAA
 TACTAGAAAGAAAATCATTGTG**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-Mlul
 ACCN: NM_173129
 Insert Size: 936 bp



OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_173129.1</u> , <u>NP_775152.1</u>
RefSeq Size:	1241 bp
RefSeq ORF:	936 bp
Locus ID:	286913
Cytogenetics:	15p12
Gene Summary:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]