

Product datasheet for **RN214517**

Adck1 (NM_001108985) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Adck1 (NM_001108985) Rat Untagged Clone
Tag: Tag Free
Symbol: Adck1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN214517 representing NM_001108985
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCCAGAAAGGCTATCAAGCTTGCTTCGTGGACCAGCATGGCTCTTGCTGCCTCCGGGGTCTACCTCT
ACAGCAACAGCTACTTGGACCCTAGTGACTTTGGCGTGT CAGGGTGGCAGAGCTGTTGCTACGACAGC
TGTCATCAGCTATGACTACCTCACCTCCCTGAGGAGTGCCCATATGGCTCAGAGGAGTATTTGCAGCGT
CGATCCCAGATCCATGATTTGTTCTGAGCTTCGATGACACCCCTCTCGGGCAGCCTCCCTGGCCAGG
TCCACAAGGCAGTATTGCATGATGGCCGACAGTGGCAGTGAAAGTCCAGCACCCAAAAGTTCAGGCTCA
AAGCTCCAAGGACATTCTCCTGATGGAGGTGCTTGTCTGGTGTGAAGCAACTTTCCAGATTTTGAA
TTCATGTGGTGGTGGATGAAGCCAAGAAGAACCTGCCTCTCGAGTTGGACTTCTAAACGAAGGGAGGA
ACGCTGAGAAAGTGGCCACATGCTCAAGCACTTTGACTTCTAAAGGTCCCCAGATCCACTGGGAGCT
GTCTACCAAGAGGGTGTCTCCTGATGGAATTTGTGGAGGGAGGTCAGGTCAACGACAGGGCCTACATGGAG
AAGAACC GGATCGATGTGAATGAGATCTCCTGCCACCTAGGCAAGATGTACAGTGAGATGATCTTTGTCA
ATGGCTTCGTACACTGTGACCCCAACCAGGCAATGTGCTGGTGGGAAGCGTCTGACACGGGCAAGGC
TGAGATTGTCTCTGGATCATGGGCTCTACCAGGTGCTCACGGAGGAGTTCCGCCTGGACTACTGCAAT
CTATGGCAGTCTCTGATCTGGACCGACTTGGACAGAGTGAAATATTACAGCCAGCGCCTGGGGGCTGCAG
ACCTCTACCCACTGTTTGCTGTATGCTGACCGCCGGTCTTGGGCTCAGTCAAACAGGGCATTGGGCA
AGCTCCAGTCTCTGCTACTGAGGACTCAGAGATTTCGCAATAATGCAGCCTGCTACCTGCCTGAGATCAGC
CAGCTCCTTAACCATGTGCCTCGCCAGATGCTGCTCATCTGAAGACCAACGATCTACTCCGAGCATCG
AGACCACCTGGGCACCGGGTCCAGCGCAGTTCTTCTCAACATGTCTCGGTGTTGTATCCGGGGCGT
GGCTGAGCATCGGAAGAGGGATGCCGGCTATTCTT CAGAAGGACCCAGATATCTTT CAGTGAGGCTTC
AGCCTGTGGCAGATTAGCCTTCATGAATCTGCTTCGAGTGAAGGCCTTGAGGCTAGCTTGTGGGTCT
CAGCTCTCTGGGCTGGCTGACTCGGGCTCCATAC**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001108985
Insert Size:	1368 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001108985.1</u> , <u>NP_001102455.1</u>
RefSeq Size:	2080 bp
RefSeq ORF:	1368 bp
Locus ID:	366698
Cytogenetics:	6q31