

## Product datasheet for RN214594

### Gipc2 (NM\_001037210) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Gipc2 (NM\_001037210) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Gipc2  
**Synonyms:** MGC116274  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Fully Sequenced ORF:** >RN214594 representing NM\_001037210  
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGCCGCTTGGGCTGCGTGGGAAGAAGAAGGCTGCCAAGTCCAAGGAGACCGCGAGGCTGGTGGAGGGCG  
 AGAGGTCCGGTGGCAGCCAGGGTGTCCCCGACCCCACTCCGGCACGGAGGCTGGTGTCCACGCACA  
 GCTGGCGCACGGCAGCGCCACGGGCCGAGTGGAGGACTTCTCCAGCATCAGCGAGCTCTACGCCAAGATC  
 GCAGGGGTGTTTGTAGATCGCGCCGTCGAGATCCTATTTGTACTTTAAACACACCTAAAATTGACATGG  
 GAAAACCTCTGGGTGGCAGCTAGTCTTGAAGATTTTCATATTTGCCACGTGAAAGGCATAAAAAAGGA  
 AGTGAACGTTTATAAATCAGAGGATTCACTTGGACTCACCATCACGGATAACGGCGTTGGCTATGCTTTT  
 ATCAAGAGAATTAAGGATGGCAGCACCATTGACTCAGTTAAAACGATCTGTGTGGGGATCATATCGAAT  
 GTATAAACGGAGAAAATATTGTTGGGTGGCGTCACTTTGAAGTCGCTAAGAAGTTAAAGGAATTGAAAA  
 GGAGGAACCTTTACTTTGCAGTTAATAGAACCAAGAAGGCATTTGAAATAGGACCAAGGTCAAAGGCC  
 GGAAAGACATCAACGGAAAAATTTGGTACTTCAAGAGGGACTCTTCGCCTGAGATCGAAAGGCCCGCTA  
 CGGTAGAAGATCTGCCCTCTGAAGCCAAAGCAAAGCGATCGAGAAGGTTGACGACCTCCTTGAAGTGA  
 CATGGGGATTGAGATATCGACTTAGCTACCACCATGTTTGAAGCTGGAAAAGACAAGAGCAACCCAGAT  
 GAGTTTGTGTCGCGCTCGACGAGACCCGGGAGACTTTGCCTTTCCTGATGAGTTCATGTTTCGACGTAT  
 GGGGAGCCATCAGTGATGTGAAACAAGGCAGAT**GA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001037210



[View online »](#)

<b>Insert Size:</b>	945 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001037210.1</a></u> , <u><a href="#">NP_001032287.1</a></u>
<b>RefSeq Size:</b>	1376 bp
<b>RefSeq ORF:</b>	945 bp
<b>Locus ID:</b>	365960
<b>UniProt ID:</b>	<u><a href="#">Q498D9</a></u>
<b>Cytogenetics:</b>	2q45