

## Product datasheet for **RN214599**

### MgII (NM\_138502) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** MgII (NM\_138502) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** MgII  
**Synonyms:** MAGL  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN214599 representing NM\_138502  
**Red**=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGCCTGAGGCAAGTTCACCCAGGCGAACTCCACAGAACGTCCTACCAGGACCTTCCTCACCTGGTCA  
ATGCGGATGGACAGTACCTCTTTGTAGTACTGGAAGCCAGTGGCACACCCAAGGCCCTCATCTTCGT  
GTCCCATGGAGCTGGGGAACACTGTGGCGTTATGACGAGCTGGCTCAGATGTTGAAGAGGCTGGACATG  
CTGGTGTGGCCATGACCATGTTGGCCATGGGCAGAGCGAGGGAGAGAGGATGGTGGTATCGGACTTCC  
AAGTTTTGTGAGAGATTTGTTGCAGCAGTGAACACCGTCCAGAAGGACTACCCCGAGGTCCCCGTCTT  
CCTCCTGGGCCACTCCATGGGCGGTGCCATCTCCATCCTAGCAGCTGCAGAGAGACCAACCCACTTTTCT  
GGCATGATCCTAATTTACCTCTGATCCTTGCCAAATCCGGAATCTGCATCGACTTTGAAGGTCTTGCTG  
CCAAACTGCTCAATTTGTCTGCCAAACATATCCTTGGGGCGCATCGACTCCAGCGTGTGTCTCGAA  
CAAGTCGGAGGTGACCTGTACAACCTCCGACCCACTCATCTGCCACGCAGGGGTGAAGGTATGCTTTGGC  
ATCCAGCTGCTGAACGCTGTCTCGAGGGTGGAGCGAGCAATGCCAGGCTGACACTGCCGTTCTGCTGC  
TGCAGGGTCTGCTGACCGCTTTGCGACAGCAAAGGTGCCTACCTGCTCATGGAATCATCCCGAGTCA  
GGACAAAACACTCAAGATGTATGAAGGTGCCTACCATGTTCTCCACAAGGAGCTTCCGGAAGTGACCAAC  
TCTGTCTCCATGAAATAAATACGTGGGTGCTCACAGGATAGCAGTGGCAGGAGCTAGGTGTCTACCCCT  
GA

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_138502  
**Insert Size:** 912 bp



[View online »](#)

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_138502.2](#), [NP\\_612511.1](#)

**RefSeq Size:** 1234 bp

**RefSeq ORF:** 912 bp

**Locus ID:** 29254

**UniProt ID:** [Q8R431](#)

**Cytogenetics:** 4q34

**Gene Summary:** serine hydrolase that converts monoglycerides to fatty acid and glycerol; may play a role in retrograde signaling at central synapses and other forms of short-range neuronal communication [RGD, Feb 2006]