

## Product datasheet for **RN214607**

### **C2cd4c (NM\_001109285) Rat Untagged Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** C2cd4c (NM\_001109285) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** C2cd4c  
**Synonyms:** Fam148c; RGD1561354  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN214607 representing NM\_001109285  
**Red**=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGAGAAAACCAACATGTGGTCTTGAACGACTTCGGGGTCTGGGAGAACGGAGCCAGCCGGGG  
TGGGTGGCAGGCTGGGACAAAACCTCCAAGGCCCTTTGTATAGCAACGTGCTGACTCCGACAAGT  
CCCCGACTTCTCATCCCGCCAAAGCTGCCCTCAGGTCCCCAGAGGCTGAGGGACAGCCGACCTGGG  
CCCTCCACCTCAGAGCAAACCTGGCTTCACTGGGCCCGTCGAGCACCTCGCAGCCCTCGGCTGCCG  
TCAAGTTGGCCTCAGAGAGCAGGAACCTGCTGAAGGCTGCCACACGACACGTGATTCAGATAGAGAGTGC  
GGAGGACTGGACAGCAGAGGAAGCCACCAACGCTGACCCAGGCACAGGGTGCCATGTCATTGCCCTCT  
GTGCCAAGGCTCAGACATCCTATGGCTTTGCCACACTGGCCGAGAGCCCCACACCAGGCGCAAGGAGT  
CCCTTTCCACAGCGAGCACGGGGCCCTGGCTCAGGTCGGATCCCCTGGTGTAGCCGCCGGCAGCAGG  
GGCTAAGGCCAATGGAGGTGATGGGGGTCCCGGAGGTTGGGGTGGCCTGATGAGCCCTAGCCGCTAC  
TTTAGTGGTGGAGAAAGTGACACAGGTTCTCTGCTGAGTCGTCGCCCTTCGGTCCCCTCTGCTGTCTC  
GCTCTGTATCGTTGCTCAAGGGCTTTGCCAGGACAGCCAGGCCAAAGTGAGCCAGCTGAAACAGTCGGT  
AGGTCCGCATGGCTCCCTGTCTGCGGATGACAGTACACCGGACACCCAGCCAGGTGTCCGTCCCGCTTG  
ACCCGCAGGGCTACCCCGAGCCAGGCCGAGTCTGGCCAGGCACCACGTGGAGAGCACACAGTACGGA  
TGGGCACCAGGGCAGCGTGAGGCTGCTGGCCGAGTATGAGGCTGCTCAGGCCGCTGCGTGTGCCCT  
GCTGGCGCCGAGGGCTGTATGACCGCCCTGTGATGCGCGGAGCATCAATTGCTGCGTGGGGCTGTGC  
CTGGTGCCGGGAAGCTGCAGAAGCAGCGGAGCACCATTATCAAAAACAGCCGCATCCCATCTTCAACG  
AGGACTTCTTCTTGACGGCCTGGGGCCAGCCAGTGTGCGTAAGCTGGCCCTCAGAATCAAGGTGGTCAA  
TAAGGTTAGCAGCCTCAAGCGGGACACACTCTTAGGGGAGGAGGAAGTCCACTGACCTCCTTGTGCC  
TTTCT**GTGA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



<b>Restriction Sites:</b>	Sgfl-Mlul
<b>ACCN:</b>	NM_001109285
<b>Insert Size:</b>	1269 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001109285.1</a></u> , <u><a href="#">NP_001102755.1</a></u>
<b>RefSeq Size:</b>	1597 bp
<b>RefSeq ORF:</b>	1269 bp
<b>Locus ID:</b>	500798
<b>Cytogenetics:</b>	7q11