

Product datasheet for RN214614

Olr1254 (NM_001001085) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Olr1254 (NM_001001085) Rat Untagged Clone
Tag: Tag Free
Symbol: Olr1254
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >RN214614 representing NM_001001085
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGACTCAGTAAATATCTCTTTGGTGAATTCATTCTGGTAGGATTAACAGACCAGCCTTATTTAC
 AAATACCATTGTTTCATTGTGTTCTAGCAATGTACCTCATCTGCATTGGGAAATGTATCTTTAATAAT
 ACTAAGTGTGCTGAATTCACCTTCACACCCCTATGACTTTTTCTCTTTAACTTGCCTTTGTAGAC
 CTCTGCTATTCTCTGTGTTCACTCCCCAAATGCTGATGAATTCATAAAAAGGAAAAATACAATTTCTT
 ACATGGAATGTATGGTCCAATCTATTTCTCTGTTTTTTGTTATTTCTGAGTGTATGTATTGACTTC
 AATGGCCTATGATCGCTATGTGGCCATCTGTAACCCACTGTTGTATAAACTTGTATGTCCTCTAAATTA
 TGTCTGAACCTAATGCTTGGTTCGTAATTCATTGCATTTCTGAGTCGGTAGCTCACACTGCATGCATGC
 TGAGATTGACCTTCTGTGATGCTAACACCATCACTACTACTTCTGTGACATTCCCCCTCTGCTTCAGCT
 TTCTTGCACGACCACACGTCATGAGCTTGTCAATTTCTGTTGTTGGGAGCATCAATATCATTGTTCCC
 ATTTCAACTATATTTATCTCCTATGGTTTTATTCTCTCCAGCATCTCCGCATGAGCTCTCTGAAGGCC
 GATCCAAAGCTTTCAGCACCTGTAGCTCACACATCATTGCTGTTTTCTGTTCTTTGGCTCAGGTGCGAT
 CAGGTATTTCAAACCTCTCAGCTGGGCTATGGATGAAGGAAAAATCTCCTGTCTTTACACCAAT
 GTGATTCCATGATAAATCCTCTACTCTATAGTTTGGGAACAAAGACATTAAGTTGCCCTGAGGAAAA
 CCCTGAAGAAGAGGAACTTTGA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-Mlul
ACCN: NM_001001085
Insert Size: 933 bp



OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001001085.1</u> , <u>NP_001001085.1</u>
RefSeq Size:	933 bp
RefSeq ORF:	933 bp
Locus ID:	405976
Cytogenetics:	8q22
Gene Summary:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]