

Product datasheet for RN214648

Smr3b (NM_012684) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Smr3b (NM_012684) Rat Untagged Clone
Tag: Tag Free
Symbol: Smr3b
Synonyms: Arp; P2-VA1; RATSMR1A; Smr1; Smr1g; Vcsa1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN214648 representing NM_012684
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGAAGTCACTGTATTTGATCTTTGGCCTGTGGATCCTTCTAGCATGCTCCAGTCAGGTGAGGGTGTCA
 GAGGCCCAAGAAGACAACATAATCCTAGAAGACAACAAGATCCTTCAACTCTTCCTCATTATCTTGGTCT
 TCAGCCTGATCCCAATGGTGGACAAATAGGAGTAAACAATCACTATACCCTTAAATCTTCAACCACCTCGT
 GTTCTTGTAAATCTTCCCGTTTTATCACTGGACCACCATTGGTTGTACAAGGTACCACTGAATATCAAT
 ATCAGTGGCAGCTAACTGCTCCAGACCCTACACCTCTAAGCAATCCTCCTACTCAACTTCATTCCACAGA
 ACAAGCAAATACAAAACAGATGCCAAAATCTCCAACACTACTGCGACTACCCAAAATTCCTACTGATATT
 TTTGAAGGTGGTGGCAAA**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-MluI

ACCN: NM_012684

Insert Size: 441 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_012684.4</u> , <u>NP_036816.4</u>
RefSeq Size:	646 bp
RefSeq ORF:	441 bp
Locus ID:	24867
UniProt ID:	<u>P13432</u>
Cytogenetics:	14p22
Gene Summary:	a prohormone-like protein [RGD, Feb 2006]