

## Product datasheet for **RN214671**

### **Ptk6 (NM\_001108968) Rat Untagged Clone**

#### Product data:

|                      |  |
|----------------------|--|
| Product Type:        | Expression Plasmids  |
| Product Name:        | Ptk6 (NM_001108968) Rat Untagged Clone   |
| Tag:                 | Tag Free   |
| Symbol:              | Ptk6   |
| Vector:              | pCMV6-Entry (PS100001)   |
| E. coli Selection:   | Kanamycin (25 ug/mL)   |
| Cell Selection:      | Neomycin   |
| Fully Sequenced ORF: | >RN214671 representing NM_001108968<br>Red=Cloning site Blue=ORF Orange=Stop codon |

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGTGTCTTTGGACAAGGCTCACCTGGGCTCTAAGTATGTGGCCCTCTGGGACTCAAAGCACGAACGG  
ACGAGGAGCTGAGCTTTACGGCAGGAGACCTCCTCCATGTTACCAGGAAGGAGGAACAGTGGTGGGC  
CACCTGCTGGATGCGGAAGGCAAGGCCTTGGCTGAGGGCTACGTGCCCCACAACCTACCTGGCTGAGAAG  
GAAACTGTGGAGTCTGAACCGTGGTTCTTTGGTTGCATCTCCCGCTCAGAGGCCATGCACAGGCTGCAGT  
CTGAGGACTACCCGAAGGTACCTTCTGATCAGAGTCAGCCAGAAGCCAGGAGCAGACTATGTCCTCTC  
CGTGTGGGATGCTGAGGCTGTTTCGACATTACAGGATCTGGAGAACAGTGCAGGAGGCTGCACCTGAAT  
GAGGTGGTATCCTTCTCCAGTCTGTCTGAGCTTGTGAACACCACAAGACCCACAACCTGTCCCCTGGCC  
TACAGCTCTCCATGGCCTGCTGGAAGCTCAAACTGAGCCCTGCCCCACTGGGATGACTGGGAGAGGCC  
GAGGGAGGAGTTCACGCTCTGTAAGAAGCTGGGGTCTGGCTACTTTGGGAGGTCTTTGAAGGGCTCTGG  
AAAGCCCTGGTCCGTGTGGCTGTGAAGGTGATCTCTAGAGACAATCTCCTGCACCAGCACACCTTCCAGG  
CTGAGATTCAAGCCATGAAGAAGCTGCGACACAAGCACATCCTGTCGCTGTACGCTGTGGCAACTGCAGG  
GGACCCAGTCTATATCATCACGGAGCTCATGCCAAGGGGAGCCTGCTGCAGCTATTGCGGGACTCGGAT  
GAGAAAGACCTGCCTATTTTGGAACTGGTGGACTTTGCATCACAGGTTGCCGAGGGCATGTGCTACCTGG  
AATCTCAGAATTACATTCACCGTACCTGGCTGCAAGGAACGTTCTCGTTACTGAGAACAATCTCTGCAA  
AGTGGGGGACTTTGGGCTGGCCAGGCTCGTCAAGGAGGACATCTACCTTTCCCGTGACCACAACGTCCTCC  
TACAAATGGACAGCACCTGAGGCACTTTCCCGAGGGCATTACTCCATCAAATCGGATGTCTGGTCTTTTG  
GTGTTCTCCTTCAATAATTTTTCAGCAGGGGCGAGATGCCCTACCCAGGCATGTCCAATCACGAAACCTT  
CCTGAGGGTGGATGCCGGCTACCGCATGCCCTGCCCTTGGAGTGCCACCCAACATACACAAGCTGATG  
CTTAGCTGTTGGAGCAGAGACCCTAAGCAGAGACCCTGCTTCAAAGACCTGTGTGAGAACTCTCAGGTA  
TCACTAGGTATGAGAACCTGGTCTGA

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online >](#)

|                               |   |
|-------------------------------|---|
| <b>Restriction Sites:</b>     | Sgfl-Mlul   |
| <b>ACCN:</b>                  | NM_001108968  |
| <b>Insert Size:</b>           | 1356 bp   |
| <b>OTI Disclaimer:</b>        | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).  |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).  |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol> |
| <b>RefSeq:</b>                | <u><a href="#">NM_001108968.1</a></u> , <u><a href="#">NP_001102438.1</a></u>   |
| <b>RefSeq Size:</b>           | 2142 bp   |
| <b>RefSeq ORF:</b>            | 1356 bp   |
| <b>Locus ID:</b>              | 366275  |
| <b>Cytogenetics:</b>          | 3q43  |