

Product datasheet for **RN214717**

Khk (NM_031855) Rat Untagged Clone

Product data:

| | |
|----------------------|---|
| Product Type: | Expression Plasmids |
| Product Name: | Khk (NM_031855) Rat Untagged Clone |
| Tag: | Tag Free |
| Symbol: | Khk |
| Synonyms: | KETHPRO |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| Cell Selection: | Neomycin |
| Fully Sequenced ORF: | >RN214717 representing NM_031855 Red=Cloning site Blue=ORF Orange=Stop codon |

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGC**C

ATGGAAGAGAAGCAGATCCTGTGCGTGGGGCTGGTGGTCTGGACATCATCAATGTGGTGGACAAATACC
CAGAGGAAGACACGGATCGCAGGTGCCTATCCCAGAGATGGCAGCGTGGAGGCAACCGTCCAACCTCTG
CACTGTGCTTTCTTGTCTGGAGCCCGCTGTGCCTTCATGGGCTCGCTGGCCCATGGCCATGTTGCCGAC
TTCTTGGTGGCCGACTTCAGGCCGAGGGGTGTGGATGTGTCTCAAGTGGCCTGGCAGAGCCAGGGAGATA
CCCCTTGTCTCTGCTGCATCGTCAACAACCTCCAATGGCTCCCCTACCATATTCTCTACGACACGAACT
GCCAGATGTGTCTGCTAAGGACTTTGAGAAGGTCGATCTGACCCGGTTCAAGTGGATCCATTGAGGGC
CGGAATGCATCGGAACAGGTAAGATGCTACAGCGGATAGAACAGTACAATGCCACGCAGCCTCTGCAGC
AGAAGGTCCGGGTGTCCGTGGAGATAGAGAAGCCCCGAGAGGAACTCTCCAGCTGTTCCGGCTATGGAGA
GGTGGTGTGTCAGCAAAGATGTGGCCAAGCACCTGGGGTCCGGTCAGCAGGGGAGGCCCTGAAGGGC
TTGTACAGTCGTGTGAAGAAAGGGGCTACGCTCATCTGTGCCTGGCTGAGGAGGGAGCCGATGCCCTGG
GCCCGACGGCCAGCTGCTCCACTCAGATGCCTCCACCACCCCGAGTAGTAGACTCTCGGGGCTGG
AGACACCTTCAATGCCTCTGCATCTTCAGCCTCTCCAAGGAAACAGCATGCAGGAGGCCCTGAGATTC
GGGTGCCAGGTGGCTGGCAAGAAGTGTGGCTTGACAGGGGTTTGATGGCATTGTGTGA

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA

| | |
|--------------------|------------|
| Restriction Sites: | Sgfl-RsrII |
| ACCN: | NM_031855 |
| Insert Size: | 897 bp |



[View online »](#)

| | |
|-------------------------------|---|
| OTI Disclaimer: | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP). |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | <u>NM_031855.3</u> , <u>NP_114061.1</u> |
| RefSeq Size: | 1307 bp |
| RefSeq ORF: | 897 bp |
| Locus ID: | 25659 |
| UniProt ID: | <u>Q02974</u> |
| Cytogenetics: | 6q14 |
| Gene Summary: | an enzyme that catalyses the phosphorylation of fructose to fructose-1-phosphate [RGD, Feb 2006] |