

Product datasheet for RN214725

Olr1413 (NM_214821) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
 Product Name: Olr1413 (NM_214821) Rat Untagged Clone
 Tag: Tag Free
 Symbol: Olr1413
 Synonyms: Olfr30; TPCR05
 Mammalian Cell Selection: Neomycin
 Vector: pCMV6-Entry (PS100001)
 E. coli Selection: Kanamycin (25 ug/mL)
 Fully Sequenced ORF: >RN214725 representing NM_214821
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGAGCACCGAGCCCTGGAACCACTCCATCACTGGCTTTGTGCTGACGAGCCTTTTTGACAACAGCCAGA
 CACACCTGTTCTGTTTAGCATGGTGGTGTGGTCTACATCCTTGCTATGGTGGGCAACACAGCCATGGT
 CCTTTTGATTTGGATGGACCCCGACTCCACCCCCATGACTTCCTCCTCAGCCAGTTGTCCTTTCTG
 GACATCTTCTTTACCTCAGTCACTGTCCCAAGATGATGGTAGGCTTTCTTTGGTTGGACAAGCATCT
 CTTTCGGAGGCTGTGGGCACAGATGTTTTCTTCATGTTCTCGGGGCTGCAGAGTGCCTCCTGCTGGC
 CCTCATGGCTTACGATCGTTATGTGGCTATCTGCAACCCTCTGCGCTACCCAGTGCCTCATGAGCGGCCGG
 GTGTGCCCTGCTCATGGTCGTGGCCTCCTGGTTGGGAGGATCCCTCAACGCCTCCATTGACTTCTCTGA
 CCCTTCAGTTCCCCTACTGTGGATCACGGAAGATCTCCCACTTCTTGTGAGGTGCCCTCGCTGCTGAT
 GCTGGCCTGTGCAGACACTGAAGCCTATGAGCAGGTAATTTGTGACAGGCGTGGTGGTCTCCTGGTG
 CCCATTACATTACTGCTCTTATGCCCTCATCCTGGCTGCTGTGCTCCGAATGCACTCTGCGGAGG
 GGCGTCAGAAGGCCCTAGCCACATGCTCCTCTCACCTGACAGTCGTAATCTCTTCTATGGGCCCTTGT
 CTACACCTACATGTTACCTGCTTCTATCACTCACCAGGCAAGACGACGTAGTATCCGCTTTTTACACC
 GTTCTCACACCCATGCTTAACCTGTCTACAGCCTCAGGAACAAGGAAGTTACAGGGGCAATGAAGA
 AAGCAGTGGGGAGGTGTAGAGTCAGTAGGAATGCC**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
 ACCN: NM_214821



[View online »](#)

Insert Size:	948 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_214821.1</u> , <u>NP_999986.1</u>
RefSeq Size:	948 bp
RefSeq ORF:	948 bp
Locus ID:	252887
Cytogenetics:	10q22
Gene Summary:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]