

## Product datasheet for **RN214776**

### **Rab12 (NM\_013017) Rat Untagged Clone**

#### Product data:

Product Type: Expression Plasmids  
Product Name: Rab12 (NM\_013017) Rat Untagged Clone  
Tag: Tag Free  
Symbol: Rab12  
Mammalian Cell Selection: Neomycin  
Vector: pCMV6-Entry (PS100001)  
E. coli Selection: Kanamycin (25 ug/mL)  
Fully Sequenced ORF: >RN214776 representing NM\_013017

Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGGATCCGAGCGCCGCGCTGCACCGCGGCCCGCGGGCGGGCCTGGGCGGGTCTCCCCGGCGCTGT  
CCGGCGGCCAGGCCCGTTCGGAGGAAGCAGCCCCCAGGCCGGCCGACTTCAAGCTGCAGGTATCATCAT  
CGGCTCCCGCGGGCTGGGCAAGACCAGCCTGATGGAGCGCTTACCGACGACACGTTCTGCGAGGCCCTGC  
AAGTCCACCGTGGGTGTTGACTTTAAAAATCAAGACTGTAGAACTAAGAGGGAAGAAAATTAGATTACAGA  
TCTGGGACACAGCAGGGCAGGAGCGATTCAACAGCATTACCTCAGCCTATTACAGAAGTGCCAAGGGGAT  
CATACTCGTGTATGACATCACCAAGAAAGAGACGTTTCGATGACTTGCCAAAGTGGATGAAGATGATCGAT  
AAGTATGCTTCAGAAGACGCTGAGCTTCTCCTGGTCGAAATAAGCTGGACTGTAAACAGACAGAGAAA  
TCTCCAGGCAGCAAGGAGAAAAGTTTGACAGCAGATAACCGGGATGCGGTTCTGCGAAGCCAGTGCCAA  
GGACAATTTCAACGTGGATGAGATCTTTCTGAAACTCGTTGATGATATTTTGAAAAAGATGCCTCTCGAT  
GTTTTGAGAAGCGAGTTATCCAATAGCATCCTCTCTACAACCAGAGCCTGAGATCCCACCAGAGCTGC  
CTCCACCCAGACCACACGTCCGATGCTGT**TGA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI  
ACCN: NM\_013017  
Insert Size: 732 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_013017.1</a> , <a href="#">NP_037149.1</a>
<b>RefSeq Size:</b>	1984 bp
<b>RefSeq ORF:</b>	732 bp
<b>Locus ID:</b>	25530
<b>UniProt ID:</b>	<a href="#">P35284</a>
<b>Cytogenetics:</b>	9q37
<b>Gene Summary:</b>	may play a role in vesicle-mediated transport; may regulate neurotransmitter release [RGD, Feb 2006]