

## Product datasheet for **RN214836**

### Stkld1 (NM\_001107822) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Stkld1 (NM\_001107822) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Stkld1  
**Synonyms:** RGD1307355  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN214836 representing NM\_001107822  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGCAGAAATTTCTCCCATCGACCAGAGGTCCAACCTCAGAGCCCTTAACAAGCTTCTGACAATGCCAGAGG  
ATCAGCTAGGCCCTGCCGTGGCCACAGAGCTGGTAGAGGAGGTCATCACCATCATAAAGCAGCATGAGAA  
GGTCCTGGATATCTGCTCAGCGCTGCTCCCTGCTGCTGCGTGTCTGGCCAAGCACTGGCGCATGAC  
TCAGAAGCTGAGATCCCGAAGAGCAGCTTGATCACCTCCTTCCTGATGAATGCCTTGAGAGCCACCCTA  
ACTCTGAAGGCCTTGTCTTGTGACCTACAATGTGCTCAGCATCATCTCTAGCCAAGGACAGATCTCAGA  
GGAGCTAGAAGAGAAGAGATTGTTTGTGCTTGGCCAAAAGCACCTGGAGCGCTTCCAAGAGGACAGGGAC  
ATCTGCCCTCGGTATCCTGAGCCTGCTCTGGTCCCTCCTGGTAGATGTTGACACTGTGAACAAAGAGCCCT  
TGGAGAACCTCTCAGTCTCAGTCACCTGGGTGCTGGCTTTTCATCAAGAGGACATGGAAGTAACAGAGGC  
TGGCTGCGCAGTCTCTGGCTTCTGTCTTGTGCTGATGGGCTGCATAAAGGAGGTCAGTTTGAGCAGGTGGT  
GTGCTGCTGCTGAGAAGCATCCAGCTATGCTCCTGGCAGAGTACTTCTGGTGAACAATGCACTCCGCGGT  
TGGCCAGCCTCGAAAGGTGTCTGAACCTGGTGGCCTTCCGAACAGTAGTACTGGAAGAGGGCAGCAGTGG  
CCTCCACCTCGTCCAAGACATCTACCAGCTCTACAAGGATGACCCCTGAGGTGGTGGAGAACATCTGCATG  
CTGTTGGCCAACTGGCCTCCTACAAGGAGATCCTGCCAGAGATGATGTCTGGAGGCATCAAGACCTAG  
TCCAGGTGATCCGGGGCGCTTTACCTCCAGTCTGGAGCTGATTTCTTATGCTGATGAGATACTTCAGGC  
ATTGGAAGCAGCTACACATCCTAGCCCCAGGAGGATCAGCTTGAACCTCCTGCAATGCTGGAAGCCTCC  
AAGGTCTCCCACTGCAGGCAGATCCTCTCCTCAGGCC**TGA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001107822



|                               |   |
|-------------------------------|---|
| <b>Insert Size:</b>           | 1092 bp   |
| <b>OTI Disclaimer:</b>        | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).  |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).  |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol> |
| <b>RefSeq:</b>                | <u>NM_001107822.1, NP_001101292.1</u>   |
| <b>RefSeq Size:</b>           | 2270 bp   |
| <b>RefSeq ORF:</b>            | 1092 bp   |
| <b>Locus ID:</b>              | 311825  |
| <b>Cytogenetics:</b>          | 3p13  |