

Product datasheet for **RN214869**

Stk19 (NM_001013197) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Stk19 (NM_001013197) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Stk19
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
Fully Sequenced ORF:	>RN214869 representing NM_001013197 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGAATCGTAAGAGACGTCTCTGGCGTCTGAGGTTTTGGAGTGAAGCGACGGCGGTGCCGGGGCCCC
TGAGAGCGGATCCTCTAAGGGCGGAAGCAGGTTTCGGCACGCGAGGCCATCGAAGAGCTCGTGCAGCTGTT
CCCACGGGGCTGTTTCGAGGACGCGCTGCCGCCATCGCGCTGAGGACTCAAGTGTACAGTCTGGTGCCG
GACAGGACGGTGGCCGACCTGCAGCTGAAGGAGCTTCAAGATCGGGGAGAGATCCAAGTCATCCAGCTAG
GCTTTGACTTGGATGCCACGGAATTGTCTTCACTGAAGACTACAGGACCAGGTCCTCAAGGCCTGCGA
TGGCCGACCATGCGCTGGGGCGGTGCAGAAGTTCCTGGCCTCGGTGCTCCAGCCTGTGGGGACCTTAGT
TTCCAGCAGGACCAGATGACACAGACTTACGGCTTCAGGGACTCGGAGATCACGCACCTGGTGAATGCTG
GAGTCTCACTGTCCGTGATGCTGGAAGCTGGTGGCTGGCTGTGCCCGAGCTGGGAGATTCATTAAGTG
CTTTGTTAAAGGACGCCAGGCTGTACTGAGCATGGTGCAGGCAAGTACCGGGAGCTTCCCCTGTGC
GAGCTCTGGGCCGAGAGTCCCTTCAGCAGTGGCTAGGTCTGCCTACCATGTGCATGACCTCATTG
GGCACAGCTGGTGGACTGTGTCCACCACCTTCTGGAACCTCCTTCGCTGCCAGATACGTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001013197
Insert Size:	765 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001013197.2</u> , <u>NP_001013215.2</u>
RefSeq Size:	1161 bp
RefSeq ORF:	765 bp
Locus ID:	361800
Cytogenetics:	20p12