

Product datasheet for **RN214891**

Rbbp9 (NM_019219) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Rbbp9 (NM_019219) Rat Untagged Clone
Tag: Tag Free
Symbol: Rbbp9
Synonyms: MGC124617
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN214891 representing NM_019219
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGCGTCCCCTACCAAGGCAGTGATTGTTCTGGGAACGGAGGGCGGGATGTGGCCACCCACGGCTGGT
ACGGCTGGGTGAGAAAGGGCTGGAGCAGATTCCTGGTTCCAGTGTGGCTAAAAACATGCCTGACCC
AATTACCGCTCGAGAGAGCATCTGGCTGCCCTTCATGGAGACAGAAGTGCCTGTGATGAGAAGCCATC
ATCATAGGCCACAGTTCGGGGCCATCGCAGCCATGAGGTATGCAGAGACACATCAGGTATACGCTCTCA
TATTGGTGTCTGCATACACATCAGACTTGGGAGATGAAAATGAGCGTGCAAGTGGGTACTTCAGCCGCC
CTGGCAGTGGGAGAAGATCAAGGCCAACTGCCCTCACATTATACAGTTTGGCTCTACTGATGACCCCTTC
CTTCCATGGAAGGAACAACAAGAAGTGGCAGATAGGCTGGACGCCAACTGTACAAATTCCTGACCGTG
GTCACCTTCAGAACACAGAGTTCATGAACTGATTAGAGTGGTGAAGTCTATGCTGACTCCTGCTGTGA
A

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_019219
Insert Size: 561 bp
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_019219.2</u> , <u>NP_062092.2</u>
RefSeq Size:	1898 bp
RefSeq ORF:	561 bp
Locus ID:	29459
UniProt ID:	<u>O88350</u>
Cytogenetics:	3q41
Gene Summary:	binds retinoblastoma (Rb); may play a role in transformation [RGD, Feb 2006]