

Product datasheet for **RN214896**

Mapk1ip1 (NM_001122782) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Mapk1ip1 (NM_001122782) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Mapk1ip1
Synonyms:	RGD1565431
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN214896 representing NM_001122782 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGTATCCCTTCCTCCCTCCGGCTAGACTGTTACCAGGATCTCCAGCACCCCTTTCTTCTTCTGGACCAT
CATGTCCCAACCTAGTGGTCCTTATCCAGGCCCTGCTGTACGAGTCCCTGGCCACACAAAGCCATATGT
TACAACAAATATGCCCTTCCAGAGCTACCTAGGCCAAATAGTGACACACAGATCCAGCTGGTCCTTTA
GGTACACAGGGATCCATGTCTTCTGGACCTTGGCGCCAGGAATGGGAGGGCAGCATCCCAATGTGCCAT
ATCTATTTCCAGAGCCATCTCCCACTCCTCCCCTGCCAGTGTGAGGAGCACCACCTGTTGCATGGGTCAC
AGTGCCACCAGGAGCCTGGAAACCACCGGCATCCTATCCTTCTCCTGAAGCTTCATATCCACCCAGGA
CTCCAGCCTTCTCCGAATAATTCTACCAACTGCCTCCAGGACCTTCTGGTGCTTCACCAATGCCTGGTA
ACCTTCATAAAGTGAATGAAATGCCAGGTGGCTCCCATTTCTGATTCGTCACCCGAGAGCACCTTGAT
GACCACTGGGCAAAAGAAGCAGCTGAAGGGAGACAACAAAGCCATCAAGAGGAGACGTTACAAGAAGAAA
AGCAAGCGAGTAACCTGGGGAGACATCAAGACTCTAACTCACAAGCTGAGACATTGGGGAAACAACAAG
GACACAATACTACTGACCCCAAAATGATGCTGCTTTGTTAATGACTATGCTGCATGTTAATTCTCAGCA
GGAACACGAAGGGTCCAA**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001122782
Insert Size:	792 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001122782.1, NP_001116254.1</u>
RefSeq Size:	1445 bp
RefSeq ORF:	792 bp
Locus ID:	499280
Cytogenetics:	1q41