

## Product datasheet for **RN214920**

### Irf3 (NM\_001006969) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Irf3 (NM\_001006969) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Irf3  
**Synonyms:** MGC94729  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN214920 representing NM\_001006969  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGGAACCCGAAACCGCTGATTTTGCCTTGGCTGGTGTCTCAGCTGGACCTGGGACAGCTGAAAGGTG  
TGGCCTGGCTGGACGAGAGCCGTACAAAGTTCCAGATCCCATGGAAGCATGGCTTACGACAGGACGCACA  
GATGGCTGACTTTGGCATCTTCCAGGCCTGGGCGGAAGCCAGTGGTGCCTACACCCAGGGAAGGATAAG  
CCAGACCTGTCAACCTGGAAGAGGAATTTCCGGTCAGCCCTGAACCGGAAAGAAGTGTTCGATTAGCTG  
AGGACCGGAGCAAGGACCCTTTTACCCTCATAAAGTGTATGAGTTTGTGACTCCAGGAGGCAAGGGA  
CTTTGTACATCTGGACACCTCCCCTGACACCAATGGCAAAGCAGTCTGTCTGATCACCAGGAAGACCTC  
TTGGAATTACTCGATCACATGGCCTTGGGACCCCTCCCAGATGAGGGGTCTCAGACCTGCCTATTGCTT  
CTGATCCTTCTCAACCACCACTAAGCCCCATTGTAACAACCTCCCAAACCCAGCACCCAGGAAAAACCC  
ACTAAGGCAGCTGCTAGCTGAGGAACAGTGGGAGTTCGAGGTGACTGCCTTCTACCGAGGCCGCCAGGTC  
TTCCAGCAGACACTCTTTTCCCCGGGGCCCTGCGGCTGGTGGCAGCAGCTTGACAACGGGACACTGC  
CCTGGCAGCCAGTCACCCTGCCAGACCCCTGAGGAGTTTCTGACAGACAGGCTTGTGAGGGAGTATGTGAG  
GCAGGTAICTCAAGGGGCTGGGCAAGGGGCTGGTCTGTGGCGGGCAGGGCAGTGCCTCTGGGCCAGCGC  
CTAGGCCACTCGCATTCTTCTGGGCCCTGGGTGAGGAGCTGTTCCAGACAGTGGGAGAGGGCCTGATG  
GAGAGTCCCCAAGGACAAGAACGGAGTCGTGTTCCGACCTCAGGCCCTTTGTGGCAGATCTCATTGCCTT  
CATGGAAGGAAGCAGACATTTCCACGATAACACTCTGTGGTTCTGTGTGGGGGAATCGTGGCCCCAGGAC  
CAGCCGTGGTCAAGAGGCTTGTGATGGTCAAGGTTGTTCTACATGTCTTAAGGAGCTGTTAGAGATGG  
CCCGGGAAGGGGAGCCTCATCACTGAAAACCGTGGACTTGCACATCTCCAACAGCCAGCCGATCTCCCT  
TACCTCTGACCAGTACAAGGCCCTGCCTCCAGGACTTGGTGAAGACATGGACTTCCAGGCCACTGGAGAA  
ACCTGA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

<b>Restriction Sites:</b>	Sgfl-Mlul
<b>ACCN:</b>	NM_001006969
<b>Insert Size:</b>	1266 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001006969.1</a></u> , <u><a href="#">NP_001006970.1</a></u>
<b>RefSeq Size:</b>	1471 bp
<b>RefSeq ORF:</b>	1266 bp
<b>Locus ID:</b>	292892
<b>Cytogenetics:</b>	1q22