

## Product datasheet for **RN214962**

### **Rcan1 (NM\_153724) Rat Untagged Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Rcan1 (NM\_153724) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Rcan1  
**Synonyms:** Dscr1; Mcip1  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN214962 representing NM\_153724  
**Red**=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGCATTTTAGGGACTTTAACTACAATTTAGCTCCCTGATTGCCTGTGTGGCAAACGGTGATGTCTTCA  
GCGAAAGTGAGACCAGGGCCAAATTTGAATCCCTCTTCAGGACGTATGACAAGGACATCACCTCCAGTA  
TTTTAAGAGCTTCAAGCGTGTCCGAATAAACTTCAGCAACCCCTGTCTGCGGGCAGCCAGGCTGCAG  
CTGCATAAGACAGAGTTCTGGGAAGGAGATGAAACTGTACTTTGCCAGACTTACACATAGGAAGTT  
CACACCTGGCTCCGCCAATCCAGACAAACAGTTTCTCATCTCTCCCTGCCTCTCCACCTGTTGGCTG  
GAAACAAGTAGAAGACGCTACCCAGTCATAAACTACGATCTTTTATATGCTATCTCCAAGCTGGGACCA  
GGAGAGAAGTACGAGCTACACGCAGCGACAGACACCACGCCAGCGTGGTGGTCCACGTGTGTGAGAGTG  
ACCAAGAGGAGGAGGAGGAGGAGGAGGAGATGGAGAGAATGAAGAGACCCAAGCCAAAAATCATCCA  
GACTCGGAGGCCAGAGTACACCCATCCACCTCAG**TGA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_153724  
**Insert Size:** 600 bp  
**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_153724.2</a></u> , <u><a href="#">NP_714946.2</a></u>
<b>RefSeq Size:</b>	2165 bp
<b>RefSeq ORF:</b>	600 bp
<b>Locus ID:</b>	266766
<b>UniProt ID:</b>	<u><a href="#">Q6IN33</a></u>
<b>Cytogenetics:</b>	11q11
<b>Gene Summary:</b>	Inhibits calcineurin-dependent transcriptional responses by binding to the catalytic domain of calcineurin A. Could play a role during central nervous system development. [UniProtKB/Swiss-Prot Function]