

Product datasheet for **RN215002**

Wdfy1 (NM_001008331) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Wdfy1 (NM_001008331) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Wdfy1
Synonyms:	MGC93914
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN215002 representing NM_001008331 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGTGCACGCGCAGCGGGAACATGCTTGCCGTCACCTTCTCTCGTCATGGGCTTCCTGTCTGCAGTACG
 ATCTGGATACTCAGCATGCATTTGTGGCGATTACTCGGGACAGATCACCTGCTGAAGCTGGAGCAGAA
 CTCCTGCTCAGTTATCACCACTCAAGGGCCATGAAGGTAGTATCGCCTGCCTTTGGTGGGACCTATC
 CAGCGTTACTCTCTCAGGAGCATCCGACAACAGCATCATCATGTGGGACATCGGGGTCGAAAGGCC
 GAACACTCCTCCTTCAGGGCCACCATGACAGGGTGCAGTCACTGTGCTACCTGCAGCTCACCAGGCAGCT
 CGTCTCCTGCTCGCGGATGGCGGAATTGCGGTGTGGAACATGGATGTTAGCCGAGAAGAGGCTCCTCAG
 TGGTTAGAAAGCGATTCTGTGAGAAGTGTGAGCAGCCCTTCTTCTGGAATATAAAGCAGATGTGGGACA
 CGAAGACACTGGGGCTGAGGCAGGAAGTGTGGGCAGGCCGTCTGTGGGAAGTGCAGCAGCAAGAGGTCCA
 GCTACCTGTCTATGGGCT**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001008331
Insert Size:	582 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001008331.1, NP_001008332.1</u>
RefSeq Size:	1441 bp
RefSeq ORF:	582 bp
Locus ID:	301549
Cytogenetics:	9q34