

Product datasheet for RN215067

Rpia (NM_001108632) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
 Product Name: Rpia (NM_001108632) Rat Untagged Clone
 Tag: Tag Free
 Symbol: Rpia
 Mammalian Cell Selection: Neomycin
 Vector: pCMV6-Entry (PS100001)
 E. coli Selection: Kanamycin (25 ug/mL)
 Fully Sequenced ORF: >RN215067 representing NM_001108632
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCAGCGCCCCGGGCCCTTCAGCACCCCTACGGGCGGGTCTGGCCCCACTGCCCGGGCGGGCCGGG
 GCGCGGCTCGGGTGGGGCGGGAACTGGGCCTCTCGGGCTCCACGTGCAGTTGCCGGGGCGTAC
 ACAGTCGGAGACCCGCGGCGACAAGGGAGGCTCGAGCGGGCGGGCCCGCCCTCCACGATGTCCAAG
 GCCGAGGAGGCTAAGAACTGGCGAGCCACACGGCAGTGGAGAACCACGTGAAGAATAACCAAGTGTGG
 GAATTGGAAGTGGTTCTACAATTGTCCATGCTGTGCAGCGAATAGCTGAAAGAGTGAAGCAAGAGAACCT
 GGACCTCATCTGCATCCCCACATCTTCCAGGCCAGGCAGCTTATCCTACAGTATGGCTTAACCCTCAGT
 GACCTGGATCAACACCCGGAGATCGACCTTGCCATCGATGGTGTGATGAAGTGGATGCTGAGCTCAATC
 TCATCAAGGGTGGTGGAGGCTGCCTGACCCAGGAGAAGATTGTGGCTGGTTATGCCAGTCGCTTCATTGT
 GATCGCTGATTCAGGAAAGATTCAAAGAACCTTGGGGATCGGTGGCACAAGGGGATTCTATTGAAGTC
 ATCCCAATGGCCTATGTCCAGTGAGTCGAGCTGTGACCCAGAAGTTGGGGGTGAGGTGGAACCTCGAA
 TGGCTGTCAACAAGGCAGGGCCTGTGGTGACAGATAATGGGAATTTATCCTGGACTGGAAGTTCGATCG
 GGTACATAAATGGGGTGAAGTGAACACAGCTATCAAGATGACTCCAGGCGTGGTGGACACAGGACTTTTC
 ATCAACATGGCAGAGAGAGTCTACTTCGGGATGCAGGATGGCTCAGTGAACGTGAGGGAGAAGCCCTTCT
 AA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-MluI
 ACCN: NM_001108632
 Insert Size: 912 bp



[View online »](#)

| | |
|-------------------------------|---|
| OTI Disclaimer: | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP). |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | <u>NM_001108632.1</u> , <u>NP_001102102.1</u> |
| RefSeq Size: | 1812 bp |
| RefSeq ORF: | 912 bp |
| Locus ID: | 362383 |
| Cytogenetics: | 4q31 |