

## Product datasheet for **RN215205**

### Rpp25 (NM\_001012124) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Rpp25 (NM\_001012124) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Rpp25  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN215205 representing NM\_001012124  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGAGAACTTCCGTAAGGTGCGCTCGGAGGAGGCGCCGGGGGACGGTGATGAGGGTGGCAGCCCGA  
GCTCTGGCCCTTTTCAGATCTGGCGCTGGTGTACACATGCGGGTCAAGGAGGGCAGCAAGATCCG  
CAACTGCTGGCCTTCGCCACCGCTAGCATGGCGCAGCCAGCCAGCGCGCCATCGTCTTCAGCGGCTGC  
GGCCGGGGCCACCACCAAGACCGTCACGTGCGCCGAGATTCTCAAGCGCCGCTGGCGGGCTTACACCAGG  
TCACGCGGCTGCGCTACAGGAGCGTGCGGAGGTGTGGCAGAGCCTCCCGCTGGGCCACACCGGGTCA  
GACGCCAGCGATCCGGCCGCGTCTCAGTGTACTTAAGAATGTACCGAGCCTCGCCATCTACTTTCA  
AAGGACGCACTGGATCCACGACAGCTTGGCTACCAGCCCCGAACCTCAGTCCTGGTCTTCGTCCCTC  
CTACCGTGTGACGTCGAAGAGGAGCCTGGGGGAATCTGCAGCTGGAGAAGGCACAGCTAAGCGGTCTCA  
GCCTGAGCCAGAGGCTGAGAATGAGGACCGGACTGCC**TGA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001012124  
**Insert Size:** 600 bp  
**OTI Disclaimer:**

Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

|                               |   |
|-------------------------------|---|
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).  |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol> |
| <b>RefSeq:</b>                | <u><a href="#">NM_001012124.1</a></u> , <u><a href="#">NP_001012124.1</a></u>   |
| <b>RefSeq Size:</b>           | 1405 bp   |
| <b>RefSeq ORF:</b>            | 600 bp  |
| <b>Locus ID:</b>              | 315705  |
| <b>UniProt ID:</b>            | <u><a href="#">Q5PPN2</a></u>   |
| <b>Cytogenetics:</b>          | 8q24  |
| <b>Gene Summary:</b>          | Component of ribonuclease P, a ribonucleoprotein complex that generates mature tRNA molecules by cleaving their 5'-ends. Also a component of the MRP ribonuclease complex, which cleaves pre-rRNA sequences.[UniProtKB/Swiss-Prot Function]   |