

Product datasheet for RN215231

Tmem169 (NM_001109574) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Tmem169 (NM_001109574) Rat Untagged Clone
Tag: Tag Free
Symbol: Tmem169
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >RN215231 representing NM_001109574
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAAGACTCAGCACCGGTGGGAAGCCAGGGTCAGCTCCAAGCCCCACCACGGCTCACTGAGGAGGG
CGGTGGCTGCTGTCTGGCCCTGGATGGGAATCCACATTGGGCCGTAGGAAGAAGAGGAGGAAAGACTC
TCGCCCAGAATCCATCATCATCTATCGGTCAGATAATGAGAAGACGGATGAGGAACCTGACGAATCTGAA
GCTGGAGATCGGCCTAAGGAAGAAGAGGGGAGGACTTCCTGGACTACCTGGAGATGATGGTATGTGGA
ACATGCCCTCTGGACAGTCGCTATGTCACACTAAGTGGGACCATCACACGAGGGAAGAAAAAGGGTCAGAT
GGTGGATATCCACGTCACCTCTGACAGAGAAGGAGCTGCAGGAGCTGACCAAACCCAGGAGCTGTCAAGG
GAAGCAGCACCCGAGGGAAGAAGGGCCTGTGAGGTGGGGCAGACCAGGGGCCACATGTGGTCTCTGGA
CACTAGTCTGCCTGCCTGTGGTCTTTGTCTCTCTTTGTGTCTTTCTACTATGGCACCATTACCTG
GTACAACATCTTTCTTGTGTAATGAAGAAAGGACCTTCTGGCATAAGATTTCTGCTGCCGTGCCTC
ATTCTTCTATCCAGTGTCTCATGATGGCTATGGCCTCTCCCTGGGGCTTTATGCTGCCGTGGTCCAGC
TCTCATGGTCTGGGCAGCATGGTGGCGGGCTGCCTGTGACATGGAGAAGGGCTTCTGTGGCTGGCTCTG
CAGCAAGCTGGGTCTGGAGGACTGTTCTCCCTACAGCATTGTGGAGCTGCTGGAGTCTGATAATATCTCC
AGCAATCTCTCCAACAAGGACCCATCCCGGAAGTAGAAACCTCTACTGT**TAG**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-MluI
ACCN: NM_001109574
Insert Size: 894 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001109574.1, NP_001103044.1</u>
RefSeq Size:	2942 bp
RefSeq ORF:	894 bp
Locus ID:	690294
Cytogenetics:	9q33