

## Product datasheet for **RN215245**

### Rragc (NM\_001048184) Rat Untagged Clone

#### Product data:

Product Type: Expression Plasmids  
Product Name: Rragc (NM\_001048184) Rat Untagged Clone  
Tag: Tag Free  
Symbol: Rragc  
Vector: pCMV6-Entry (PS100001)  
E. coli Selection: Kanamycin (25 ug/mL)  
Cell Selection: Neomycin  
Fully Sequenced ORF: >RN215245 representing NM\_001048184  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGTCCCTGCAATACGGGGCGGAGGAGACGCCCTGGCCGGCAGTTACGGCGCAGCCGACTCATTCCCCA  
AGGACTTCGGCTACGGCGTGGAGGAGGAGGAAGAAGAGGCAGCGCGGGTGGCGCGGAGGCGCGGGGGC  
CGGCGCGGCTGCGGCCCGGGAGGCGGTGACAGCTCCAAGCCAAGGATCTTGCTTATGGGACTACGGCGC  
AGCGGCAAGTCTCCATTAGAAGGTGGTGTTCATAAGATGTACCCAATGAGACTCTCTTTTGGAAA  
GTACCAACAAGATTTATAAAGATGACATTTCCAACAGCTCCTTTGTCAACTTCCAGATTTGGGATTTTCC  
TGGGCAGATGGACTTCTTCGACCAACCTTTGACTATGAAATGATCTTCAGGGGGACAGGAGCATTGATA  
TATGTCATTGATGCACAGGATGACTACATGGAGGCTTTAACTCGACTTCACATCACTGTATCTAAGGCC  
ACAAAGTGAACCCAGACATGAACCTTTGAGGTTTTTATTCACAAAGTTGATGGTCTGTCTGATGATCACA  
GATAGAAACACAGAGGGACATTCATCAAAGGGCCAATGATGACCTGGCAGATGCTGGGCTGGAAAAGCTC  
CATCTTAGCTTTTATTTGACTAGCATCTATGACCACTCAATATTTGAAGCTTTCAGTAAGGTGGTACAGA  
AACTCATCCCGCAACTGCCACTTTGGAAAACCTATTAATATCTTTATATCAAATTCAGGTATTGAAA  
AGCTTTTCTCTTTGATGTTGTCAGCAAAATCTACATTGCAACAGACAGCTCCCCTGTGGATATGAGTCT  
TACGAGCTCTGCTGTGACATGATTGATGTTGTAATTGATGTGTCCTGTATATATGGATTAAGGAAGATG  
GAAGTGGAAAGTGCTTATGACAAAGAATCGATGGCCATTATCAAGCTGAACAATACAACTGTCTTTATTT  
AAAGGAAGTCACTAAATTTTGGCACTGGTTTGCATTCTTAGGGAAGAGAGCTTTGAACGAAAAGTTTTG  
ATAGACTACAACCTTCACTGTTTCCGAAAAGCTATCCATGAGGTGTTTGGGTGGGTGTGACTTCTCACA  
GAAGCTGTAGTCATCAGACCAGTCTCCAAGCCTGAAAGCCTTGGCACACAATGGCACACCTCGGAATGC  
CAT**TAG**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



[View online »](#)

<b>ACCN:</b>	NM_001048184
<b>Insert Size:</b>	1197 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_001048184.1</a></u> , <u><a href="#">NP_001041649.1</a></u>
<b>RefSeq Size:</b>	2548 bp
<b>RefSeq ORF:</b>	1197 bp
<b>Locus ID:</b>	298514
<b>Cytogenetics:</b>	5q36